

XRF Analysis of Polymers Additives, Initiator Residues and Toxic Elements

Arkady Buman and Dr. Kai Behrens
Bruker AXS Inc.

Welcome Today's Speakers



Arkady Buman
Product Manager XRF
Madison, WI, USA



Dr. Kai Behrens
Global Product Manager XRF
Karlsruhe, Germany

Webinar: XRF Analysis of Polymers Topics




- Analytical tasks in polymer and plastics production
 - Process control
 - Quality control
- Analytical technology: WDXRF
- Results
 - Additives
 - Initiator residues
 - Impurities
 - Regulated toxic elements
- Ready-to-analyze solutions
- Q&A



Polymers and Plastics

Homopolymers



- PE, PP, PVC and PS are very popular polymers and have replaced a lot of traditional materials, such as metal or wood, due to longer lifetime and weight, better shaping and coloring characteristics
- ABS, PET, and Teflon are replacing other materials, such as ceramics and glass, in technologically advanced products

Polymers and Plastics Additives




- Thermal stabilizers (processing, long term)
- UV stabilizers
- Pigments
- Carbon black
- Processing aids
- Anti-block
- Anti-statics
- Slip agents
- Flame retardants


Polymers and Plastics Stabilizer and Processing Aids




- Metal deactivators form stable complexes with metal ions thus improving stability



- Calciumstearate



- Magnesiumstearate



- Mg, P, Cu and Ca are indicators for additives
- Mg, Al, Ti, Ni, Mo are typical residues from initiators
- F, Cl, Fe, Zn are impurities from production processes
- Pb, Hg, Br, Cd, Cr are regulated substances (RoHS, WEEE), other regulations such as consumer product safety and end of life of vehicles (ELV) also regulate other elements such as Sb

S8 TIGER WDXRF Spectrometer Series




S8 TIGER⁽⁽⁽ WDXRF spectrometers are the latest models on the market with most modern technology!

Available in three basic configurations:

- S8 TIGER 1K
- S8 TIGER 3K
- S8 TIGER 4K

S8 TIGER

Sequential WDXRF Spectrometer



Analytical flexibility:

- 4 kW excitation
20 – 60 kV
- 5 – 170 mA
- 10 beam filters
- 4 collimators
- 8 crystals
- 2 detectors

WDXRF Spectrometer Series: S8 TIGER 4K




Superior analytical performance with clear benefits by design:

- Best in the market:
- max. 170mA - 4 kW for best light element analysis
- No compressed air
- Dust-sealed spectrometer cabinet
- Smallest footprint (~0.8 m²)
- Reduced helium consumption
- Reduced cooling water consumption, independent of short fluctuations

S8 TIGER

Superior Analytical Performance



Unrivalled analytical flexibility:

- 10-position beam filter changer
- 4-position collimator changer
(6 different divergence settings)
- 8-position crystal changer
(15 different crystals)

Application-optimized crystals:

- XS-CEM: long-term stable, temperature-independent
- XS-GE-C: plus 40% intensity for P, plus 20% intensity for S
- XS-PET-C: plus 20% intensity for Al
- XS-B: plus 100% intensity for B


S8 TIGER

X-ray Technology – High Performance



- Highest excitation intensity for all elements
 - Optimized beam path with highest flexibility
 - Compact beam path with lowest distance anode – sample - detector
- Best light element excitation
 - Low kV excitation voltage for reduced background
 - Full 4 kW excitation power up to 170 mA for best element sensitivities
- Excellent precision
 - Direct loading of the sample in the reference position, no inaccurate positioning in comparison to turret systems
 - Reliable analysis of liquid samples
- No heat-up and burning of samples in comparison to conventional XRF systems
 - Efficient tube head cooling

S8 TIGER
SampleCare™



S8 TIGER SampleCare™





Four Times Protection:

- Low maintenance cost and time, high instrument uptime due to unique protection during loading and unloading:
- Two contamination shields to protect tube window and goniometer
- Dust reservoir to collect sample particles and droplets

Unique protection during measurements:

- Duraberyllium™ shield for tube window protection
- Unique vacuum seal with high transmission window for goniometer protection



S8 TIGER Polymer Applications



S8 TIGER

Analysis of Additives




- Determination of additives
 - Precise evaluation of the concentration level for the optimal dosing of valuable substances
 - Selection of tracer element for organic additives, which cannot be determined directly
 - e.g. P as marker for PS 38 or PS P-EPO
- Analysis of solids as hot pressed PUK under vacuum
- Analysis of granules in sample cups under helium
- Selection of soft excitation condition to reduce the heat dissipation from tube to sample



Additives in PE


Analysis of Phosphorus with the S8 TIGER



- Separation of P K α from neighboring peaks
- Use of high intensity and best resolution crystal XS-Ge-C

Additives in PE

Analysis of Phosphorus with the S8 TIGER



Additives in PE

Analysis of Phosphorus with the S8 TIGER




Number	Standard Name	Chemical Concentration	XRF Concentration	Absolute Deviation	Relative Deviation	LLD (PPM)
14	AB-1	1.3	1.2	-0.1		0.3
10	AB-2	15.7	16	0.3	1.7	0.3
20	AB-3	30.1	29	-1.1	-3.5	0.3
17	AB-4	45.8	42.6	-3.2	-7.1	0.3
15	AB-5	46.6	46.4	-0.2	-0.36	0.3
13	AB-6	61	60.5	-0.5	-0.75	0.3
11	AB-7	74.8	77.4	2.6	3.4	0.4
16	AB-8	98.2	97.7	-0.5	-0.48	0.5
12	AB-9	110.3	110.1	-0.2	-0.16	0.5
18	AB-10	125.8	126.9	1.1	0.87	0.5
19	AB-11	160.8	160.4	-0.4	-0.23	0.6

- Direct analysis of solid samples (disks - PUK) from a hot press or as granules
- Detection limit of 0.3 ppm (3σ , 60 s)
- Optimal detection with clear separation from other elements

Additives in PE


Analysis of Calcium with the S8 TIGER



- Separation of Ca K α from neighboring peaks
- Flat background, no background determination required

Additives in PE

Analysis of Calcium with the S8 TIGER



Additives in PE

Analysis of Calcium with the S8 TIGER



Number	Standard Name	Chemical Concentration	XRF Concentration	Absolute Deviation	Relative Deviation	LLD (PPM)
1	P1	2.3	2.5	0.2		0.5
2	p2	21.6	22.8	1.2	5.4	0.5
3	p3	50	50.1	0.1	0.24	0.6
4	p4	69.4	70.2	0.8	1.2	0.7
5	p5	95.6	94	-1.6	-1.7	0.9
6	p6	121.5	118.5	-3	-2.5	1
7	p7	139.6	140.4	0.8	0.6	1.1
8	p8	165	162.9	-2.1	-1.2	1.1
9	p9	186.2	189	2.8	1.5	1.1
10	p10	210	209.8	-0.2	-0.073	1.4
11	p11	232.7	233.6	0.9	0.37	1.4

- High accuracy allows exact determination of the mixture of polymer with the additive
- Maximum deviation 2 ppm (1% relative)
- Detection limit of 1 ppm (3 σ , 60 s)
- Less thermal stress on samples reducing the measurement time
- No background measurements required

S8 TIGER

Analysis of Additives



- Additives are easily determined by WDXRF
- Detectable elements such as P used as marker for organic substances
 - Providing chemical formula allows the direct calculation of the compound concentration
- Reducing the measurement time allows the analysis of a wide element range without stressing the sample
 - Low temperature tube head with efficient cooling reduces the heat dissipation to the sample
- Analysis of solid samples or granules without complicated sample preparation




S8 TIGER

Analysis of Initiator Residues




- Analysis of initiator residues
 - Mg beside high Al concentration in Ziegler-Natta Systems
 - Cr or Ti in 3d-element based initiator system
- Early recognition of initiator losses requires precise trace element determination
- Avoidance of high initiator concentrations in order to maintain long polymer life
- Meet regulation such as packaging directives, food safety and consumer product safety

Initiator Residues in PE Analysis of Magnesium with the S8 TIGER



- Impurity from production
- Efficiency improved by reduced Mg concentration
- LLD: 0.2 ppm

Initiator Residues in PE Analysis of Chromium with the S8 TIGER



- Calibration Std. Dev. of 0.1 ppm
- Dedicated concentration range to max. 6 ppm for optimal trace detection

Initiator Residues in PE


Analysis of Chromium with the S8 TIGER



Number	Standard Chemical Name	Concentration	XRF Concentration	Absolute Deviation	Relative Deviation	LLD (PPM)
1	CR1		0.8	0.8	0	0.2
2	CR2		3.4	3.2	-0.2	0.2
3	CR5		4.5	4.8	0.3	0.2
4	CR4		4.8	4.8	0	0.2
5	CR3		6	5.9	-0.1	0.3

- Optimal accuracy of less than 0.3 ppm deviation allows early recognition of initiator losses
- Detection limit of 0.2 ppm (3σ , 60 s)
- Excellent precision
- Dedicated trace calibration enhances WDXRF performance
- Unrivalled accuracy of WDXRF in comparison to other methods such as ICP or AAS – due to little sample preparation

Initiator Residues in PE Analysis of Titanium with the S8 TIGER



- Calibration Std. Dev. of 1 ppm
- Dedicated concentration range to max. 200 ppm

Initiator Residues in PE

Analysis of Titanium with the S8 TIGER



Number	Standard Name	Chemical Concentration	XRF Concentration	Absolute Deviation	Relative Deviation	LLD (PPM)
1	TITAN003	0.1	0	-0.1		0.2
2	TITAN009	14.9	14.4	-0.5	-3.3	0.3
3	TITAN019	22.4	22.1	-0.3	-1.1	0.4
4	TITAN023	28.8	27.1	-1.7	-5.8	0.5
5	TITAN034	40.8	39.9	-0.9	-2.3	0.6
6	TITAN038	45.3	44	-1.3	-2.9	0.7
7	TITAN048	65.9	65.9	0	-0.025	1
8	TITAN052	78.3	77.2	-1.1	-1.4	1
9	TITAN061	131.3	133.1	1.8	1.4	1

- Direct analysis of solid samples (disks - PUK) from a hot press or as granules
- Detection limit of 0.3 ppm (3σ , 60 s) for low ranges or 1 ppm for higher ranges
- 1 ppm deviation for higher concentration ranges

S8 TIGER

New X-ray Tube

High-intensity X-ray tube with Cr target
Operating at 3.3 kW

Target applications:

- Petrochemistry
- Polymers

Available for

- S8 TIGER



S8 TIGER

New Cr X-ray tube




Intensity gain (Cr over Rh tube)

- For Cl - factor of 3
- For Ca - factor of 4
- For Ti - factor of 5

- LLD in polymers
 - 0.4 ppm for Cl (0.7 ppm with Rh)
 - 0.1 ppm for Ca (0.2 ppm with Rh)
 - 0.05 ppm for Ti (0.2 ppm with Rh)

- Less stress on valuable samples due to 3.3 kW max power



S8 TIGER


Analysis of Impurities



- Impurities from the production process influences polymer quality, such as iron and other elements from steel used for the piping and machineries
- Halogens decrease polymer life
- In solid samples the detection of F is possible with WDXRF in hot pressed samples

Impurities in PE

Analysis of Fluorine with the S8 TIGER



- Separation of F K α from neighboring peaks, such as Fe L
- Use of the high resolution crystal XS-55 (multilayer optic)

Impurities in PE

Analysis of Fluorine with the S8 TIGER



Number	Standard Name	Chemical Concentration	XRF Concentration	Absolute Deviation	Relative Deviation	LLD (PPM)
1	F1	0	-1	-1		11.7
2	F2	1630	1904	274	17	12
3	F3	3418	4667	1249	37	12.9
4	F4	4378	4424	46	1	12.8
5	F5	6543	6264	-279	-4	13.3
6	F6	7444	6400	-1044	-14	13.4
7	F7	10238	10409	171	2	15.2

- Direct analysis of solid samples (disks - PUK) from a hot press in vacuum
- Analysis of granules in polymer cups is not feasible due to absorption of F K α in the foil
- Detection limit of 13 ppm
- Excellent precision


S8 TIGER

Analysis of Regulated Hazardous Elements




- Detection and analysis of regulated toxic elements are possible with WDXRF down to the sub-ppm range
- EDXRF is not capable of separating elements such as Hg, Br and Pb from each other
 - Low spectral resolution

Regulated Toxic Elements Analysis of RoHS Elements with the S8 TIGER




- Optimal resolution
- WDXRF separates neighboring lines
- Clear identification of Hg, Pb, and Br

Regulated Toxic Elements Analysis of Cadmium with the S8 TIGER



- Separation of Cd from tube lines
- (Rh $K\alpha$ for the S8 TIGER)
- with LiF 200, better with LiF 220

Regulated Toxic Elements Analysis of Cadmium with the S8 TIGER



- Calibration range up to 200 ppm
- Calib. Std. Dev. of less than 4 ppm
- High linearity

Regulated Toxic Elements Analysis of RoHS Elements with the S8 TIGER



Sample	Date	Cr (%)	Br (%)	Hg (%)	Pb (%)	Cd (%)
Average (n=83)	every 10 min	0.0479	0.0346	0.0288	0.0333	0.0017
Stdev		0.0002	0.0005	0.0001	0.0001	0.0001
COV		0.5%	1.4%	0.4%	0.4%	6.0%

Sample	Date	Cr (PPM)	Br (PPM)	Cd (PPM)	Hg (PPM)	Pb (PPM)
QC_RoHS-QUANT-ABS-6mmDisk	12/10/2008 9:18	503	227	95	401	475
QC_RoHS-QUANT-ABS-6mmDisk/1	12/10/2008 9:27	504	227	95	402	473
QC_RoHS-QUANT-ABS-6mmDisk/2	12/10/2008 9:36	503	227	93	403	474
QC_RoHS-QUANT-ABS-6mmDisk/3	12/10/2008 9:45	504	227	93	400	474
QC_RoHS-QUANT-ABS-6mmDisk/4	12/10/2008 9:54	503	227	92	402	474
QC_RoHS-QUANT-ABS-6mmDisk/5	12/10/2008 10:03	504	226	94	401	474
QC_RoHS-QUANT-ABS-6mmDisk/6	12/10/2008 10:11	504	226	95	401	475
QC_RoHS-QUANT-ABS-6mmDisk/7	12/10/2008 10:20	503	227	94	403	475
QC_RoHS-QUANT-ABS-6mmDisk/8	12/10/2008 10:29	505	227	97	402	476
QC_RoHS-QUANT-ABS-6mmDisk/9	12/10/2008 10:38	504	226	96	402	475
QC_RoHS-QUANT-ABS-6mmDisk/10	12/10/2008 10:47	504	227	94	400	476
QC_RoHS-QUANT-ABS-6mmDisk/11	12/10/2008 10:56	503	227	94	400	473
QC_RoHS-QUANT-ABS-6mmDisk/12	12/10/2008 11:05	504	226	96	401	474
Average		504	227	94	401	475
Std.Dev.		1	0	1	1	1
Rel.Std.Dev.		0.14%	0.20%	1.32%	0.26%	0.17%

- Excellent precision of less than 1.5 % relative in ppm levels

S8 TIGER WDXRF Performance



- Accuracy and precision fit to demands for the analysis of
 - Additives
 - Initiator residues
 - Impurities
 - Regulated elements
- Simple sample preparation
 - For integration into process control regime and quality checks close to production
- Detection of traces to sub-ppm with standard instrument configuration
- Ready to analyze solutions available

S8 TIGER Ready-To-Analyze Solutions



GEO-QUANT T



GEO-QUANT M



GEO-QUANT F



METAL-QUANT



CEMENT-QUANT



POLYMER-QUANT A



PETRO-QUANT



RoHS-QUANT



QUANT EXPRESS

POLYMER-QUANT A for the S8 TIGER



POLYMER-QUANT is the ready-to-analyze solution for additives in polymers:

- Covering more than 8 elements
- In polymer disks and granules
- Contains certified standard materials from NIST as acceptance test sample

POLYMER-QUANT A APPLICATIONS



- Contains all relevant additive elements and initiator residues in polymers
- Covers the analysis of virgin polymers with typical concentration ranges
- In all CH-based polymer types with the variable alpha model (PE, PP, PET, ABS, ...)


A periodic table of elements with a light gray background and rounded corners. The elements are arranged in their standard periodic layout. The table is divided into several sections: the main body (rows 1-7), the lanthanide series (row 8), and the actinide series (row 9). The elements are color-coded: H and He are white; Li, Be, B, C, N, O, F, Ne are light gray; Na, Mg, Al, Si, P, S, Cl, Ar are white; K, Ca, Sc, Ti, V, Cr, Mn, Fe, Co, Ni, Cu, Zn, Ga, Ge, As, Se, Br, Kr are light gray; Rb, Sr, Y, Zr, Nb, Mo, Tc, Ru, Rh, Pd, Ag, Cd, In, Sn, Sb, Te, I, Xe are light gray; Cs, Ba, La, Hf, Ta, W, Re, Os, Ir, Pt, Au, Hg, Tl, Pb, Bi, Po, At, Rn are light gray; Fr, Ra, Ac are light gray; Ce, Pr, Nd, Pm, Sm, Eu, Gd, Tb, Dy, Ho, Er, Tm, Yb, Lu are light gray; Th, Pa, U, Np, Pu, Am, Cm, Bk, Cf, Es, Fm, Md, No are light gray. The elements Ti, Zr, and Hf are highlighted with a white background.

POLYMER-QUANT A

Covered Concentration Range



	Mg (ppm)	Al (ppm)	Si (ppm)	P (ppm)	Cl (ppm)	Ca (ppm)	Ti (ppm)	Zr (ppm)
MIN	2	26	5	1	10	0	0	3.5
MAX	495	260	1198	210	1250	125	99	76



- Covered polymer types: PE, PP, PET, ABS, ...
- Combination with QUANT-EXPRESS for PVC and similar polymers
- Combination with RoHS-QUANT ABS
- Typical measurement: 14 min (S8 TIGER 4 K)
Adjustable to requirements

POLYMER-QUANT A Standards



POLYMER-QUANT A

- Contains calibration standards (Polyethylene and Polypropylene) as granules
- NIST certified reference material as acceptance test samples
- Drift correction samples
- Operators Manual
- Installation CD

S8 TIGER Ready-To-Analyze Solutions



GEO-QUANT T



GEO-QUANT M



GEO-QUANT F



METAL-QUANT



CEMENT-QUANT



POLYMER-QUANT A



PETRO-QUANT



RoHS-QUANT



QUANT EXPRESS

RoHS-QUANT ABS S8 TIGER



- Is the analytical ready-made solution for RoHS for quick & easy qualitative screening and quantitative analysis of regulated substances in polymers of electrical and electronic devices for the
 - S2 RANGER
 - S8 TIGER



RoHS - Lead-free

Limits of elements and compounds



Elements	Defined Limits	Use
Lead	1000 ppm	paint, solder, PVC stabilizer, plastic pigment
Mercury	1000 ppm	batteries, contacts in microswitches, lightning equipments, plastic pigment
Cadmium	100 ppm	stabilizer and pigment in plastics, semiconductors, batteries
Chromium	1000 ppm (Hexavalent)	anti-corrosion coatings, plasticizers, paints, plastic processing
Bromine	1000 ppm (as PBB and PBDE in total)	fire retardent in plastics and polymers

RoHS Quantitative Screening Mandatory Detection Limits



- Basic requirements for the safe and error-free identification of elements and the quantitative analysis of regulated elements are the resolution and detection power of the XRF instrument.
- Not every XRF instrument is capable of distinguishing between Br and Hg or to detect safely Cd in RoHS samples.
- Mandatory detection limits in ppm for XRF used for RoHS compliance screening acc. to IEC/ACEA:

	Polymer	Metals	Electronics
Cd	15	30	30
Pb	30	60	60
Hg	30	60	60
Br	15	30	
Cr	30	60	30

S8 TIGER RoHS-QUANT



RoHS-QUANT is based on fast and simple sample preparation:

Benefits:

- Quick and instant results
- Flexible for different kind of samples
- Low costs for sample preparation equipment or accessories

RoHS-QUANT ABS S8 TIGER



- 10 ABS polymer standards (100 g granules each)
- 10 ABS polymer standards prepared as 2- and 6-mm disks
- 1 QC sample (100 g, 2- and 6-mm disks)
- CRM BR680 PE standard (47 g and 2-mm disk)
- CRM BR681 PE standard (47 g and 2-mm disk)

S8 TIGER Ready-To-Analyze Solutions



GEO-QUANT T



GEO-QUANT M



GEO-QUANT F



METAL-QUANT



CEMENT-QUANT



POLYMER-QUANT A



PETRO-QUANT




RoHS-QUANT



QUANT EXPRESS

Refinery Products



PETRO-QUANT

Solutions for the Petrochemical Industry



- Norm-compliant to ASTM, EN, ISO, DIN
- Sulfur and lead in automotive fuels, lubricants
- Additives and stabilizers in oils, greases, ...
- Wear metals to detect engine debris
- Impurities of initiators in oils and polymers
- Additives in polymers and plastics
- Toxic elements according to RoHS/WEEE
- Pet coke / coal anodes
- Alternative fuels in cement production
- ...



PETRO-QUANT

Unique Solutions for Petrochemicals



H																	He	
Li	Be											B	C	N	O	F	Ne	
Na	Mg											Al	Si	P	S	Cl	Ar	
K	Ca	Sc	Ti	V	Cr	Mn	Fe	Co	Ni	Cu	Zn	Ga	Ge	As	Se	Br	Kr	
Rb	Sr	Y	Zr	Nb	Mo	Tc	Ru	Rh	Pd	Ag	Cd	In	Sn	Sb	Te	I	Xe	
Cs	Ba	La	Hf	Ta	W	Re	Os	Ir	Pt	Au	Hg	Tl	Pb	Bi	Po	At	Rn	
Fr	Ra	Ac																
			Ce	Pr	Nd	Pm	Sm	Eu	Gd	Tb	Dy	Ho	Er	Tm	Yb	Lu		
			Th	Pa	U	Np	Pu	Am	Cm	Bk	Cf	Es	Fm	Md	No			

Universal petrochemical calibration for up to 30 elements in hydrocarbon-based matrices – straight out of the box
NEW: 4 additional elements: Co, As, Tl, Bi as additives and contaminants

Bruker AXS

Ready-to-Analyze Solutions - Benefits

POLYMER-QUANT A, ROHS-QUANT and PETRO-QUANT are fully calibrated for the analysis of additives, initiator residues, impurities and regulated elements in plastic matrices.

Benefits:

- Time saving – ready to start
- Routine analysis right from the beginning
- No XRF expert needed to do calibrations
- Minimum operator training
- QC procedure integrated
- Drift correction for the life of the instrument established

Standards included, flexibility integrated

- Quick adaptation to own applications
- Upgrade path



Q&A



Any questions?

Please type any questions you may have for our speakers in the [Q&A panel](#) and click Send.



How did we do?

When you exit the webinar, please fill out our [evaluation survey](#) to let us know. We appreciate your feedback.




Arkady Buman
Product Manager XRF
Madison, WI, USA



Dr. Kai Behrens
Global Product Manager XRF
Karlsruhe, Germany

Thank you!

Past and Future XRF Webinars



The screenshot shows the Bruker website's navigation bar with links for Industries, Products, Service, News, Events, and About us. Below the navigation is a hero section with a photo of people in a meeting and the text "Globally interacting with experts". To the right of the photo is a "Related links" section with "News" and "Events" links, and a "Contact" button. Below the hero section is a breadcrumb trail: Home • Service • Education & Training • Webinars • XRF. The main content area features a section for "XRF Webinars" with the sub-heading "Live and On-Demand". A large box highlights the URL: www.bruker.com/service/education-training/webinars/xrf.html. Below this are two columns of webinar listings, each with a "Webinar" column and a "Content" column.

Industries **Products** **Service** **News** **Events** **About us**

BRUKER

Globally interacting with experts

Related links

- ▶ News
- ▶ Events

[Contact](#)

Home • Service • Education & Training • Webinars • XRF

XRF Webinars
Live and On-Demand

Available at
www.bruker.com/service/education-training/webinars/xrf.html

Webinar	Content	Webinar	Content
Nov 20, 2012	This one-hour live webinar demonstrates the capabilities of TXRF for trace element analysis of air (aerosol and filters), land (soils and sediments), and water (fresh and effluents). Learn about level of detection, ease-of-use, regulated analysis, and advantages of TXRF in various, everyday environmental applications. Register now	Dec 13, 2010	Join Mike Beauchaine and Andrew Toms as they present the capabilities of TXRF and ICP-MS for trace element analysis of wastewater. Learn about the combined advantages of the two techniques for routine, industrial and mobile lab testing. View recording Download slides
TXRF for Trace Element Analysis of Air, Land and Water		Trace Element Analysis of Industrial Wastewater and Sewage with TXRF and ICP-MS	

Like what you learned in this webinar?

Subscribe to Bruker's *FIRST Newsletter* to get webinar announcements, fascinating articles, and analytical X-ray news delivered right to your inbox!

Subscribe at:

<https://www.bruker.com/about-us/register.html>

Upcoming Events – US and Canada Visit Us!

APS – American Physical Society – Baltimore, MD – Mar 18-22

Pittcon – Philadelphia, PA – Mar 18-21

MRS Spring – San Francisco, CA - Apr 1-5


AEROMAT – Bellevue, WA - Apr 2-5

Cast Expo – St. Louis, MO - Apr 6-9

ACS Spring – New Orleans, LA - Apr 7-11

ICDD XRF Clinic – Newtown Square, PA - Apr 29 - May 3

AISTech – Pittsburgh, PA – May 6-9



Innovation with Integrity