


Forensics Analysis with X-Ray Diffraction and X-Ray Fluorescence



Bruker AXS Inc., Madison, WI, USA
April 16, 2013



Welcome




Today's Topics

- Introduction to XRD
- Examples of forensic applications with XRD
- Introduction to TXRF, learn about X-ray fluorescence analysis
- Trace elemental analysis with TXRF
- Introduction to μ XRF for forensics
- Application examples of μ XRF
- Question & Answer

Speakers



Nathan Henderson, Ph.D.
Applications Scientist – D2 PHASER
Madison, WI, USA



Michael Beauchaine
Business Development Manager Americas - XRFi
Madison, WI, USA

Forensics Analysis with X-ray Diffraction





Forensics Analysis with X-ray Diffraction



- Basic Principles of X-Ray Diffraction
- Typical Applications
- Benchtop Diffraction with the D2 PHASER
- Selected Forensics Applications
 - Pigment Analysis
 - Pharmaceuticals
 - Explosives


A Typical Powder Diffraction Experiment



What are X-rays? ...and why do we use them?



Electromagnetic spectrum



Wavelength of X-rays on same magnitude as atomic spacing
Probe to determine atomic order and crystal structure


The Interaction of X-rays and Matter

Scattering from Multiple Atoms




Constructive Interference (peak!)

Deconstructive Interference (background)




Bragg's Law



$$\sin \theta = \frac{\lambda/2}{d} = \frac{\lambda}{2d}$$

$$2d \sin \theta = \lambda$$

Bragg Brentano Geometry




Sample Holders



- Various cavity dimensions
- Low-background holders
- Back-loading holders
- Clay slides
- Air-sensitive holders



Powder X-ray Diffraction Applications



- Qualitative Phase Analysis (Phase ID)
- Quantitative Phase Analysis
- Crystallite Size & Microstrain Analysis
- Structure Determination and Refinement

Powder diffraction *CAN* tell you: Structure, Crystallinity, Phase Composition

Powder diffraction is *NOT* an elemental analysis technique

D2 PHASER

XRD in a Benchtop Design



Smart Design

- Small footprint
- Innovative goniometer
- Bragg-Brentano geometry
- Sample rotation (up to 80 rpm)
- Choice of X-ray tubes and detectors

System Mobility

- No pre-installation
- Standard wall sockets supported
- Integrated monitor
- Internal cooling system



Pigment Analysis



Red Pigments

- Iron(III) oxide, Fe_2O_3
- Lead(II,IV) tetroxide, Pb_3O_4
- Mercuric sulfide, HgS

Orange Pigments

- Lead(II) chromate/Lead(II) oxide,
 $\text{PbCrO}_4/\text{PbO}$

Yellow Pigments

- Lead(II) chromate, PbCrO_4
- Cadmium sulfide, CdS

Green Pigments

- Chromium(III) oxide hydrate,
 $\text{Cr}_2\text{O}_3 \cdot 1.5 \text{H}_2\text{O}$

Blue Pigments

- Cobalt(II) stannate, CoSnO_3
- Iron(II,III) hexacyanoferrate(II,III),
 $\text{Fe}_7(\text{CN})_{18}$

Violet Pigments

- Cobalt(II) phosphate hydrate,
 $\text{Co}_3(\text{PO}_4)_2 \cdot 8\text{H}_2\text{O}$

Black Pigments


- Iron(II,III) oxide, Fe_3O_4
- Lead sulfide, PbS
- Carbon black, C

White Pigments

- Titanium dioxide, TiO_2
- Zinc oxide, ZnO
- Lead carbonate hydroxide,
 $2\text{PbCO}_3/\text{Pb}(\text{OH})_2$

Pigment Analysis

Discrimination between TiO_2 and ZnO



Pigment Analysis


Various Forms of TiO_2



Anatase

Tetragonal


Spacegroup: $I 41/a m d (141)$



Rutile


Tetragonal

Spacegroup: $P 42/m n m (136)$



Pigment Analysis

Quantification of Minor Phases



Pharmaceuticals




- Identification of active pharmaceutical ingredients (API), excipients, adulterants, fillers, and binders
- Common Adulterants
 - Sugars
 - Sucrose
 - Lactose
 - Dextrose
 - Mannitol
 - Caffeine
 - Acetaminophen
 - Methylsulfonylmethane (MSM)
 - Dextromethorphan
 - Corn starch
 - Sodium bicarbonate



Pharmaceuticals


Headache Powder (Aspirin and Caffeine)



Mannitol used in flavored formulations rather than lactose (common filler material)

Pharmaceuticals

Various Sweeteners (Adulterants)




Explosives and Blast Residue



- Raw Explosives
 - Positive identification of explosive materials
 - Discrimination between types of explosives
 - Pyrodex
 - Triple Seven
 - Ammonium Nitrate
- Blast Residues
 - Indication of type of explosive used
 - Oxide and halide residues
 - Aluminum chloride hexahydrate


Explosives



Introduction to TXRF



Principles X-ray Fluorescence (XRF) Spectroscopy




1. An X-ray quantum hits an inner shell electron in a (sample) atom. The electron is removed leaving the atom in an excited state
2. The missing inner shell electron is replaced by an electron from an outer shell
3. The energy difference between the inner and outer shell is balanced by the emission of a photon (fluorescence radiation)

Principles of X-ray Fluorescence (XRF) Spectroscopy




Each element shows a specific line pattern in a spectrum depending on the orbitals involved

- L→K transition = $K\alpha$ line
- M→K transition = $K\beta$ line
- M→L transition = $L\alpha$ line
- N→L transition = $L\beta$ line



Traditional EDXRF



Traditional EDXRF

Samples for common XRF spectrometry (ED and WDXRF)

- Solids (cut, polished and put into suitable shape)
- Powders (as pressed pellets, fused beads or loose powders in liquid cups)
- Liquids (in liquid cups)


Necessary sample amount: from 1 g to 10 g!!



Principles of Total Reflection X-ray Fluorescence (TXRF) Spectroscopy



Total reflection X-ray fluorescence spectroscopy



Beam angle: $1^\circ / 90^\circ$

- Samples must be prepared on a reflective media
- Polished quartz glass or polyacrylic glass disc
- Dried to a thin layer, or as a thin film or microparticle

Elements Measured by the Mo PICOFOX



H																		He
Li	Be											B	C	N	O	F		Ne
Na	Mg											Al	Si	P	S	Cl		Ar
K	Ca	Sc	Ti	V	Cr	Mn	Fe	Co	Ni	Cu	Zn	Ga	Ge	As	Se	Br		Kr
Rb	Sr	Y	Zr	Nb	Mo	Tc	Ru	Rh	Pd	Ag	Cd	In	Sn	Sb	Te	I		Xe
Cs	Ba	L	Hf	Ta	W	Re	Os	Ir	Pt	Au	Hg	Tl	Pb	Bi	Po	At		Rn
Fr	Ra	A																
		L	La	Ce	Pr	Nd	Pm	Sm	Eu	Gd	Tb	Dy	Ho	Er	Tm	Yb		Lu
		A	Ac	Th	Pa	U	Np	Pu	Am	Cm	Ek	Cf	Es	Fm	Md	No		Lr

- Impossible to measure
- Difficult to measure
- Measured using K-lines
- Measured using L-lines

The Instrument S2 PICOFOX



Benchtop TXRF spectrometer S2 PICOFOX

- Metal-ceramic X-ray tube
 - Mo anode
 - Air-cooled
 - Other tubes available
- Multilayer monochromator
- XFlash[®] silicon drift detector
 - Electro-thermally cooled
 - ≤ 149 eV @ MnKa 100 kcps
- Automatic version
 - 25-sample cassette



Application Examples



Sample Preparation



Liquid Samples - Toxicology



- fill sample in micro tube
 - add internal standard
 - homogenize
- pipette on carrier

Sample Preparation

Liquid Samples - Toxicology



- dry by heat / vacuum
 - load the instrument
 - start data acquisition

S2 PICOFOX

Trace Element Analysis - Toxicology



Introduction

- Analysis of nutrition-relevant elements (Cu, Fe, Zn, Se)
- Analysis of toxic elements (As, Pb, Hg)

Task

- Analysis of whole blood and blood serum standards

Sample preparation

- Serum
 - 1:10 dilution with water (p.a. grade)
 - Addition of Ga for internal standardization
- Whole blood
 - 1:1 dilution with water (p.a. grade)
 - Addition of Ga for internal standardization



S2 PICOFOX

Trace Element Analysis - Toxicology



- Toxicological analysis of heavy metals
- Precise results without sample digestion
- Accurate quantitation from Aluminum to Uranium
- Can analyze a variety of samples like hair, blood, urine, tissues



		Serum (ClinCheck L2)				Whole blood (Seronorm L2)			
Element	Unit	TXRF	Std. dev.	Reference ¹⁾	Std. dev.	TXRF	Std. dev.	Reference ²⁾	Std. dev.
Fe	mg/l	440	7.4	435	12	2.9	0.09	1.964	0.20
Cu	µg/l	66	2.2	62	2.1	1685	43	1562	312
Zn	µg/l	501	4.9	504	6.9	2194	118	2225	334
Se	µg/l	12	0.29	12	1.0	97	18	102	26

¹⁾: Atomic Absorption Spectroscopy

²⁾: Sector-Field Inductively-Coupled Plasma Mass Spectroscopy

Sample Preparation

Solid and Powder Samples




Solids are ground to fine particle size and resuspended for direct analysis without digestion



- fill powder in mortar
 - grind carefully ($< 50 \mu\text{m}$)
 - weigh about 20-50 mg
 - transfer to tube

Sample Preparation

Solid and Powder Samples




- suspend in detergent solution
 - add standard
- homogenize
- pipette on carrier

TXRF Applications

Authentication of Pharmaceuticals



Authentication and Purity Control of Pharmaceuticals



	Concentration (mg/kg)	
	Original	Generic
P	< LLD	25.8
Cl	2.8	16.7
K	35.2	28.9
Ca	10.5	243.8
Fe	0.85	5.4
Ni	< LLD	1.3
Zn	< LLD	1.9
Sr	0.06	0.74



- Spectra comparison after direct preparation
- Quantitative analysis after internal standardization
- ➔ **Fingerprint and authentication analysis**



Sample Preparation Microparticles



Microparticles are measured semi-quantitatively and non-destructively



- dab vacuum grease on carrier
 - pick-up some particles with a (glass) rod
 - drop particles on grease

TXRF Applications

Analysis of Glass Fragments




- 5 glass fragments (1 – 2 mm and 0.5 – 1 mm thick) are embedded in a resin block
- Fragments were too large for direct analysis
- Ground samples with agate motor then transferred to disk
- Applied discriminate analysis to understand spectral differences
- Can be applied to other small samples like paint chips
- **Develop an Inorganic Spectral Library**

TXRF Applications

Analysis of Glass Fragments



Spectral Comparison




TXRF Applications

Analysis of Glass Fragments



Spectral Comparison




TXRF Applications

Analysis of Glass Fragments



Standardless Quant. – Discriminate analysis



Introduction to μ XRF



What is μ XRF?




- Micro XRF focuses X-ray radiation to small concentrated areas of a samples by use of collimation
- This allows a position sensitive examination of:
 - Elemental distributions for non-homogenous samples of small particles, inclusions or non-regular shaped samples or
 - Thickness of composition of coatings
- Position sensitive analysis is required because a lot of materials are not homogenous and are not flat.
- μ XRF is not a destructive technique so it does not require complicated sample preparation and preserves evidence.
- Distribution analysis gives information over a large area.
- Rapid X, Y, Z stage allows for fast mapping of samples.

M4 TORNADO




M4 TORNADO

- High-end μ XRF system for
 - composition analysis and
 - layer thickness measurements
- Vacuum chamber (light elements!)
- SDD detector (up to 3 simultaneously)
- 50 W fine focus X-ray tube (Rh, Mo, W, Ag)
- Spot size 25 μm – Polycapillary Optics
- High speed, motorized xyz-stage for fast line scan and mapping
- Operated with ESPRIT software with multiple analysis and display features
- Large chamber for various sized objects




M4 TORNADO

Focusing by poly-cap lens



Poly-cap lens collects large angle of tube radiation and concentrates this to the sample in small spots



Forensic Applications Where μ XRF is Used



- Analysis of unknown and non-regular shaped samples:
 - Particles of glass and lacquer to identify cars involved in accidents
 - Particles of dust, soil, rocks, etc. found on clothing, shoes etc.
 - Different types of fibers
- Analysis of pigments on documents
 - Distinguish between different pigments
 - Detect covered or erased letters by trace analysis
- Analysis of Gun Shot Residue (GSR)
- Analysis of drugs for traces to identify their origin
- Identification of fraudulent money (paper) and of gems

Forensic Applications Analysis of Paper Money



Bank notes are high tech-products, they contain a lot of features to protect against counterfeiting

many different analytical methods are required to verify all of them

μ XRF can identify the different pigments and metal strips



Single element distributions for few elements. they show:

- The watermark is a density variation of the TiO_2 -pigment
- The metallic wire is made from Fe

Forensic Applications

Examination of Paper Money



Paper money is a high-tech product with a series of safety issues such as the used pigments that can be analyzed with XRF.



Mosaic image



Single element distributions of Cr



Fe and



Multielement overview



Ti

Measurement with
1225 x 606 pixel
25 ms/pixel

Forensic Applications

Identification of Fake Coins




The numismatic value of historic coins depends on the number of minted coins



Part of a coin of 2 Reichsmark, expected to be minted 1927, a year with low coin production

First tests are performed with sulfidic acid which generates Ag_2S (black color)



Measurement on the "2" (red spectrum) and "7" (blue spectrum) show the expected composition of Ag and Cu


The spectra show slight differences due to the Ag_2S layer

Forensic Applications


Identification of Fake Coins



The numismatic value of historic coins depends on the number of minted coins



The mapping shows that the "7" is a fake. It was a "6", was removed and replaced by a 7 which enhanced the value by approx. 2 orders of magnitude - seemingly!




Forensic Applications

Identification of Fake Gems




Fakes of pearls



These spectra of pearls shows Ca, Sr (from the seawater). But for the black pearl, there is also Ag.

This can be only from a Ag-coating of the pearl to make it more valuable.

Fakes of diamonds




Zirconia can be easily mistaken for diamond. Diamond is from Carbon, it shows no XRF line, only a large background scatter.

The diamond spectrum (red) shows some Au-lines from the support of the gem.


The zirconia spectrum (blue) has significant lines of Zr, Y and also Hf.

Forensic Applications

Gun Shot Residue – Detection of Elements of Interest



Textile with bullet hole



Distribution of S, Cu and Pb close to the bullet hole with high concentration on the hole and single particles around




Pixel: 780 x 760, 40 μm step size, 5 ms per pixel or approx. 1 h total

Forensic Applications

Gun Shot Residue – Shot Distance




Spectrum of a GSR contamination




Intensity of elements of explosives on the bullet hole depends on shot distance (drops down very fast).

This allows an estimation of shot distance.



Forensic Applications

Analysis of Pigments in an Art Object



Museum of
Fine Arts Leipzig
oil painting,
Max Klinger 1872
back side,
310 x 150 mm²

The back side of this painting on wood was covered by a paper. It was assumed that this was a older painting by Klinger which was covered by paper and than the "back side" was used again to paint.


μ XRF was able to measure through the paper to detect the distribution of pigments.

Forensic Applications

Analysis of Pigments in an Art Object



Museum of
Fine Arts Leipzig
oil painting,
Max Klinger 1872
back side ,
310 x 150 mm²



μ XRF mapping,
1000 x 500 pixel,
Ca, Cr, Fe, Co,
Cu, Hg, Pb

This examination shows that there is another painting on the back side and it seems to be a natural scenery. This result forced the removal of the covering from the back side and restoration of the painting.

Q & A




Any questions?

Please type any questions you may have for our speakers in the [Q&A panel](#) and click Send.




How did we do?

When you exit the webinar, please fill out our [evaluation survey](#) to let us know. We appreciate your feedback.



Nathan Henderson, Ph.D.
Applications Scientist –
D2 PHASER
Madison, WI, USA



Michael Beauchaine
Business Development Manager
Americas - XRFi
Madison, WI, USA

Thank you!

TXRF Webinars

Live and On-Demand



Industries

Products

Service

News

Events

About us



Globally
interacting with
experts

Related links

- ▶ News
- ▶ Events

Contact

[Home](#) • [Service](#) • [Education & Training](#) • [Webinars](#) • [XRF](#)

XRF Webinars

Live and On-Demand

Available at
www.bruker.com/service/education-training/webinars/xrf.html

Webinar

Content

Nov 20, 2012

TXRF for Trace Element Analysis of Air, Land and Water

This one-hour live webinar demonstrates the capabilities of TXRF for trace element analysis of air (aerosol and filters), land (soils and sediments), and water (fresh and effluents). Learn about level of detection, ease-of-use, regulated analysis, and advantages of TXRF in various, everyday environmental applications.
[Register now](#)

Webinar

Content

Dec 13, 2010

Trace Element Analysis of Industrial Wastewater and Sewage with TXRF and ICP-MS

Join Mike Beauchaine and Andrew Toms as they present the capabilities of TXRF and ICP-MS for trace element analysis of wastewater. Learn about the combined advantages of the two techniques for routine, industrial and mobile lab testing.
[View recording](#)
[Download slides](#)

The banner features a collage of scientific images: a circuit board, a glowing orange sphere, a test tube with blue liquid, and a pile of brown grains. The 'FIRST' logo is in large, bold, black letters with white dots. To its right, a blue bar contains the text 'FRONTIERS IN RESEARCH SCIENCE & TECHNOLOGY' in white, and 'newsletter' in white on a dark blue background. The Bruker logo, a blue atom symbol with the word 'BRUKER' in black, is on the right. A small 'ISSUE' label is in the bottom left corner.

FIRST

FRONTIERS IN RESEARCH SCIENCE & TECHNOLOGY
newsletter




Like what you learned in this webinar?

Subscribe to Bruker's ***FIRST Newsletter*** to get webinar announcements, fascinating articles, and analytical X-ray news delivered right to your inbox!

Subscribe at:

<https://www.bruker.com/about-us/register.html>



Innovation with Integrity