


# Inclusion Analysis / Steel Cleanliness Determination by Spark OES



Characterization of inclusions in steel by  
OES Pulse Discrimination Analysis (OES-PDA)




# Inclusion Analysis / Steel Cleanliness Determination by Spark OES



## Topics

- Demand & Goals
- Reference Methods for Inclusion Analysis
- Rapid Method for Inclusions & Oxygen
- Hardware - Instrument Features
- Single Spark Evaluation
- Determination of Inclusions
- Principles of the PDA/MCI Method
- Comparison of Methods
- Results of PDA/MCI
- Conclusion



## Speakers



**Georg Schick**

VP, Industrial Sales & Marketing  
Bruker Elemental  
Billerica, MA, USA



**Martin Tilleman**

Product Manager, Stationary OES  
Bruker Elemental  
Germany

# Demands & Goals

Increasing demand for higher cleanliness in steels by customers in the metal processing industries.

- Production of high purity steels
  - Process monitoring of production steps
  - Final inspection of outgoing material
- Optimization of the production process
  - Improving Slagging Practices
  - Fast determination of oxide and sulfide inclusions during the process
  - Calculation of oxygen content
- Processing of steel for high demands
  - Fast cleanliness check of incoming billets and outgoing components




**→ Determination of parameters that define Steel Cleanliness**

# Reference Method for Inclusion Analysis: SEM/EDS with Bruker Quantax 400 EDS



Scanning electron microscope with energy dispersive X-ray spectroscopy

- Universal method: differentiation of carbides, oxides, nitrides, sulfides
- Large observation area
- Imaging method
- Highest accuracy
- Surface method, low penetration depth ( $\sim 1 \mu\text{m}$ )
- Costly, long measurement time ( $\sim 3-10 \text{ h}$ )
- Highly educated operating staff




# Reference Method for Oxygen Analysis: melt extraction with G8 GALILEO



Melt extraction with carrier gas method for the determination of oxygen

- Accurate analysis of total oxygen
- Fast measurement (~80 s)
- High analysed sample mass (~1000 mg)
- Demanding sample preparation
- Limited to oxygen only



# Rapid Method for Inclusions & Oxygen: OES-PDA = MCI = Metal Cleanliness Inspection




## Inclusion characterization & oxygen determination by Optical Emission Spectrometry with Pulse Discrimination Analysis

- Complete elemental analysis
- Determination of various oxide and sulfide inclusions
- Calculation of total oxygen
- Simple sample preparation (grinding w/ SiC paper or milling)
- Fast measurement (~5 s/burn. multiple burns recommended. e.g. 5x)
- User-friendly software for "normal" OES operator
- Feasibility study advisable

A screenshot of the Bruker OES software interface. It displays a table of analysis results with columns for elements and their concentrations. The table is titled "Analyseergebnis" and contains several rows of data. Below the table, there are various fields and buttons for further analysis and reporting.

Element	Concentration	Unit
C	0.0001	%
Si	0.0001	%
Mn	0.0001	%
P	0.0001	%
S	0.0001	%
Fe	0.0001	%
Al	0.0001	%
Ni	0.0001	%
Cr	0.0001	%
Cu	0.0001	%
Zn	0.0001	%
Co	0.0001	%
Nb	0.0001	%
Mo	0.0001	%
Ti	0.0001	%
V	0.0001	%
B	0.0001	%
As	0.0001	%
Sb	0.0001	%
Bi	0.0001	%
Pb	0.0001	%
Sn	0.0001	%
Zr	0.0001	%
Hf	0.0001	%
Y	0.0001	%
La	0.0001	%
Ce	0.0001	%
Pr	0.0001	%
Nd	0.0001	%
Sm	0.0001	%
Eu	0.0001	%
Gd	0.0001	%
Tb	0.0001	%
Dy	0.0001	%
Ho	0.0001	%
Er	0.0001	%
Tm	0.0001	%
Yb	0.0001	%
Lu	0.0001	%
Sc	0.0001	%
Yt	0.0001	%
Er	0.0001	%
Tm	0.0001	%
Yb	0.0001	%
Lu	0.0001	%





# Hardware Instrument Features




## Q8 MAGELLAN Spectrometer


- Single Spark Detection (SSD) with visualization
- Unlimited parameters
  - any frequency up to 1 kHz
  - all channels available for SSD
  - free source curve design
- Highest spectral sensitivity. lowest detection limits
- Lowest dark current, Channel-Photo-Multipliers (CPM)



# Single Spark Evaluation Identification of Coincidences




Example for single spark signals  
with the Q8 MAGELLAN





# Determination of Inclusions





## Application field:

- Low and high alloyed steels with aluminium desoxidation

## Workflow calibration & measurement

- Development of inclusion method with definition of spark parameters
- Calibration of method
- Sparking samples / generating data files (CSV format)
- Import data into MCI software
- Execute MCI software
- Print out (PDF) and export (CSV) and transfer results to other databases

## Calculation steps

- Exclusion of sparks in the instable starting phase
- Detect and remove unenergetic sparks
- SD-Calculation of normalized Gaussian distribution
- Selection of oxide outliers (e.g.  $>3$  SD)
- Calculation of element concentrations
- Identification of Ca as CaS. removing from oxide examination
- Identification of different oxides (correlation of outliers, e.g. Al to Ca)
- Stoichiometric calculation of oxide and sulfide compounds
- Calculation of total oxygen
- Classification of inclusions and grain size distribution

# Comparison of Methods



	SEM/EDS	ON/(H)	OES-MCI
Capital investment (approx. k€)	550	60	80
Operating costs	High	Medium	Low
Reference method / norm compliance	Yes	Partly	No
Penetration depth (of sample), approx.	1-3 $\mu\text{m}$	Complete	10 $\mu\text{m}$
Tested area (of sample), approx.	200 $\text{mm}^2$	Complete	7 $\text{mm}^2$ )*
PDA/MCI-Measurement time, approx.	10 h	80 s )*	5 s )*
Ease-of-use (instrument)	Complex	Medium	Easy
Sample preparation	Medium	Complex	Easy
Analytical performance / value	High	Limited	Medium


)\* values given for single reading, multiple measurements necessary

# Comparison of Methods



## Validation of MCI-method by comparison with Gas Analysis

Application: low and mild alloyed steel




# Comparison of Methods



## Validation of MCI-method by comparison with SEM/EDS

Application: low and mild alloyed steel




# Comparison of Methods



## Assessment of oxide cleanliness degree

Example of calibration based on customer samples (acc. EN 10247-KOX)



# Results of PDA Method



## Reproducibility - Example ball bearing steel

Sample	Analyse No	O-total (ppm)	Al-total (ppm)	Al-insoluble (ppm)	app. K1 value (EN 10247)	PDA-value
72 - BBS $O_{REF} = 7 \text{ ppm}$	1	5.9	88	5.5	1.9	35100
	2	6.4	88	6.1	2.0	37527
	3	5.4	86	5.1	1.7	31491
	4	6.3	87	6.1	1.9	36122
	5	7.9	89	7.6	2.5	45566
	<b>Mean</b>	<b>6.4</b>	<b>88</b>	<b>6.1</b>	<b>2.0</b>	<b>37161</b>
	<b>SD or RSD (PDA)</b>	<b>0.9</b>	<b>1.1</b>	<b>0.9</b>	<b>0.3</b>	<b>14</b>
N14 - BBS $O_{REF} = 12 \text{ ppm}$	1	13.1	54	13.3	4.2	77534
	2	13.7	53	13.8	4.2	78343
	3	10.3	50	10.2	3.2	60312
	4	12.9	55	12.8	4.0	73774
	5	13.5	57	13.4	4.2	77671
	<b>Mean</b>	<b>12.7</b>	<b>54</b>	<b>12.7</b>	<b>4.0</b>	<b>73527</b>
	<b>SD or RSD (PDA)</b>	<b>1.4</b>	<b>2.6</b>	<b>1.4</b>	<b>0.4</b>	<b>10</b>



# Results of PDA Method



## Reproducibility - Example reference material, low alloy steel

Sample	Analyse No	O-total (ppm)	Al-total (ppm)	Al-insoluble (ppm)	app. K1 value (EN 10247)	PDA-value
RM - V1 *) O <sub>REF</sub> = 19 ppm	1	26.8	312	26.4	8.8	163033
	2	17.7	311	15.3	5.7	106909
	3	44.7	364	44.0	14.7	273353
	4	46.0	349	45.3	15.1	280605
	5	15.1	310	14.6	4.9	91548
	<b>Mean</b>	<b>30.1</b>	<b>329</b>	<b>29.1</b>	<b>9.8</b>	<b>183090</b>
	<b>SD or RSD (PDA)</b>	<b>14.6</b>	<b>25.5</b>	<b>14.9</b>	<b>4.8</b>	<b>49</b>
RM - V3 O <sub>REF</sub> = 17 ppm	1	19.4	1059	21.1	6.1	113989
	2	17.9	1062	19.6	5.6	103371
	3	19.3	1043	19.3	6.0	111988
	4	19.8	1061	21.3	6.2	114449
	5	17.6	1052	19.1	5.5	102769
	<b>Mean</b>	<b>18.8</b>	<b>1055</b>	<b>20.1</b>	<b>5.9</b>	<b>109313</b>
	<b>SD or RSD (PDA)</b>	<b>1.0</b>	<b>8.0</b>	<b>1.0</b>	<b>0.3</b>	<b>5</b>

\*) Inhomogeneous inclusion distribution!

# Results of PDA Method Presentation in the MCI Software



## Software Screenshot

-	22,4	115,0	27,0	260,2	13,6	111,1	18,8	-	161,0	-	-	-	-	-	-	-	-	-	-
11,9	34,2	18,1	-	862,0	115,4	248,3	21,9	469,1	-	56,6	-	-	-	-	-	-	-	-	-
8,9	-	19,3	39,2	502,3	101,8	116,8	28,2	220,6	-	-	-	-	-	-	-	-	-	-	-
11,9	3,1	71,6	19,9	461,6	-	268,2	-	506,6	100,3	-	-	-	-	-	-	-	-	-	-
18,9	18,0	14,4	-	398,2	135,8	121,4	40,7	229,2	-	-	-	-	-	-	-	-	-	-	-
-	-	59,6	-	-	-	31,5	-	-	83,4	-	-	-	-	-	-	-	-	-	-
-	9,3	188,4	-	-	-	240,2	31,3	-	-	717,5	-	-	-	-	-	-	-	-	-
-	-	25,9	9,2	-	-	-	43,8	-	36,2	-	-	-	-	-	-	-	-	-	-
-	-	50,6	0,7	-	-	-	65,7	-	70,8	-	-	-	-	-	-	-	-	-	-
-	8,1	33,7	-	-	-	-	12,5	-	47,2	-	-	-	-	-	-	-	-	-	-
5,9	-	95,1	-	-	-	3,4	15,6	-	133,1	-	-	-	-	-	-	-	-	-	-
-	-	26,5	-	-	-	13,5	43,8	-	37,1	-	-	-	-	-	-	-	-	-	-
1,0	1,9	40,3	-	-	54,3	-	21,9	-	56,5	-	-	-	-	-	-	-	-	-	-
-	13,0	101,1	-	-	-	62,1	72,0	-	141,6	-	-	-	-	-	-	-	-	-	-
-	3,7	68,6	-	-	-	-	78,2	-	96,1	-	-	-	-	-	-	-	-	-	-
3,0	-	32,5	-	-	-	-	34,4	-	45,5	-	-	-	-	-	-	-	-	-	-
-	17,4	29,5	17,0	79,2	122,2	6,7	-	-	41,3	-	-	-	-	-	-	-	-	-	-
-	11,2	-	-	-	237,6	-	-	-	-	-	-	-	-	-	-	-	-	-	-
2,0	11,8	89,1	20,6	72,4	74,7	111,1	75,1	-	124,7	-	-	-	-	-	-	-	-	-	-
2,0	-	27,7	11,3	303,2	47,5	166,1	-	313,8	38,8	-	-	-	-	-	-	-	-	-	-

	Ti	Mg	Ca	Zr	Mn	Si	Al	S	AlO	CaO	MgO	CaAlO	MgAlO	MgCaAlO	ZrO	TiO	SiO	MnO	
Statistics																			
Average	1112	1637	1705	14629	970	644													
Standard deviation	116	276	695	311	2000	177	2756												
Minimum	649	1090	1434	1700	14598	962	4853	641											
Maximum	101	246	219	305	1971	165	485	128											
Standard deviation	25	49	404	18	17	34	205	12	173	383	41	21	8	0	18	25	31	14	
Concentrations (ppm)																			
Intensity	26141	93318	117024	17009	34269	51551	1522084	32674											
Standard deviation	0,3	0,6	8,7	0,3	6,2	3,5	21,0	1,4	35,9	9,8	0,8	3,9	0,8		0,4	0,5	6,8	6,6	
Minimum		0,4	3,5				18,7												
Maximum	18,8	5,8	35,0	11,3	2415,8	206,2	406,4	31											
Standard deviation																			
Approximate total																			
Size																			
1										161	21				18	24			
2									167	220	20	2	8						
3									10			19					31	14	
4																			
5																			

Ignored sparks (Fe):	0	PDA code total:	132615,7	PDA code > 1µm:	124849,3
Ignored sparks incl. for PDA-O:	19	Clearness degree:	5,9	No. of inclusions in spark area > 1µm:	512

- Statistics
- Grain size distribution
- Total oxygen

# Results of PDA Method Presentation in the MCI Software



PDF print-out



Analysis of a ball bearing steel (Standard quality)

Q8 MAGELLAN									BRUKER		15.04.2013				
<b>Sample name</b>	<b>B140-1</b>														
Method name	Fe Inc														
No. of ignored sparks (Fe)	0														
No. of ignored CaS incl. for PDA-O	44														
PDA code total	42881,3														
PDA code >1 µm	38564,6														
No. of inclusions in spark area > 1µm	383														
No. of inclusion (per 1000mm²)	2272,7														
Cleanliness degree	1,2														
PDA Oxygen (ppm)	11,9														
<b>Statistics</b>	<b>Ti</b>	<b>Mg</b>	<b>Ca</b>	<b>Zr</b>	<b>Mn</b>	<b>Si</b>	<b>Al</b>	<b>S</b>							
Average	745	919	926	1719	25751	5785	1185	852							
Standard deviation	670	291	225	266	3041	594	722	353							
Median	691	885	897	1702	25642	5764	1059	732							
Median std. deviation	87	207	136	248	2774	563	116	116							
Outlier	168	84	61	35	20	34	261	435							
<b>Concentrations (ppm)</b>															
Outlier	4,9	1,2	1,2	0,5	24,8	10,4	12,0	35,5							
Oxygen		0,8	0,5				10,6								
Approximate total	26,6	2,5	3,8	12,2	4428,9	2758,3	45,5	94,3							
<b>Inclusion analysis</b>	<b>AIO</b>	<b>CaO</b>	<b>MgO</b>	<b>CaAlO</b>	<b>MgAlO</b>	<b>MgCaAlO</b>	<b>ZrO</b>	<b>TiO</b>	<b>SiO</b>	<b>MnO</b>	<b>TiAlO</b>	<b>MnSiO</b>	<b>CaSAIO</b>	<b>CaS</b>	<b>MnS</b>
No.	174	51	51	10	33	0	35	148	32	18	20	2	24	20	14
Outlier	13,1	0,7	1,0	0,4	5,2		0,7	3,7	20,7	15,7	5,9	2,7	6,9	2,5	14,8
<b>Grain size</b>															
Class 1	31	42	33				35	108							
Class 2	124	9	18	9	23			38		14			6	13	
Class 3	19			1	10			2	32	18	5	2	18	7	14
Class 4											1				
Class 5															

# Results of PDA Method Presentation in the MCI Software



## Classification & Presentation in Ternary Diagrams



# Conclusion



- OES-MCI reduces effort and costs, exceptionally in comparison to other methods
- Quick complementary method for material inspection
- Allows process control
- Final product assessment for steelworks and processing plants
- Easy extension of typical OES element/bulk analysis

# New Application Report

## *Steel Cleanliness Analysis by OES-PDA*



### Application Report

## Steel Cleanliness Analysis by OES-PDA

Rapid Quality Assurance Check with  
Q8 Magellan and MCI-Software

#### Abstract

The present production of steels shows demanding requirements for quality, especially the needs for homogeneity and cleanliness are increasing. The key for premium steel products are best cleanliness degrees and lowest levels of Oxygen.

The classic way to determine non-metallic inclusions is usually done with metallographic microscope and scanning electron microscope with EDX analysis (SEM-EDX). These methods are slow and cost-intensive applications. Using optical emission spectroscopy with pulse discrimination analysis (OES-PDA) means fast and efficient determination of the cleanliness degree and Oxygen level. During the melting process these quality indicators are available and allow the process control in real-time. In addition a fast quality check of the final product is one of the main advantages of this method.

For this purpose Bruker Elemental has developed the software package MCI, "Metal-Cleanliness-Inspection". The MCI-Software evaluates on statistical basis the sparks, which are caused by nonmetallic inclusions.

The main features are the determination of:

- Total oxygen content,
- Oxide and sulfide inclusions,
- Grain size distribution of different oxide types (e.g. alumina, Ca-aluminate, Mg-spinel)
- Cleanliness degree in correlation to reference methods (e.g. EN 10247, ASTM E45)

# Bruker Elemental - Billerica Demo Lab



k  
O  
E  
S  
&  
G  
A  
S




Visit us at AISTech 2013, booth #2035  
Pittsburgh, PA, USA on May 6-9



THE IRON & STEEL TECHNOLOGY CONFERENCE AND EXPOSITION 6-9 MAY 2013 • DAVID L. LAWRENCE CONVENTION CENTER • PITTSBURGH, PA, USA



[Home](#) | [Conference](#) | [Exposition](#) | [Events at a Glance](#) | [Housing & Travel](#) | [Registration](#) | [Contact](#)





# Q & A



## Any questions?

Please type any questions you may have for our speakers in the [Q&A panel](#) and click Send.



**Georg Schick**  
VP, Industrial Sales & Marketing  
Bruker Elemental  
Billerica, MA. USA

## How did we do?

When you exit the webinar, please fill out our [evaluation survey](#) to let us know. We appreciate your feedback.



**Martin Tilleman**  
Product Manager, Stationary OES  
Bruker Elemental  
Germany

**Thank you!**

# OES and CS/ONH Webinars Live and On-Demand



[Industries](#) [Products](#) [Service](#) [News](#) [Events](#) [About us](#)



Globally  
interacting with  
experts

#### Related links

- [News](#)
- [Events](#)

[Contact](#)

[Home](#) • [Service](#) • [Education & Training](#) • [Webinars](#) • [OES and CS/ONH](#)

## OES and CS/ONH Analysis Webinars

Live and On-Demand

Available at  
[www.bruker.com/service/education-training/webinars/oes-and-csonh.html](http://www.bruker.com/service/education-training/webinars/oes-and-csonh.html)

Webinar	Content
Apr 23, 2013	Join Bruker Elemental panelists Georg Schick and Martin Tilleman for this one-hour webinar in which they highlight Spark OES as a simple and easy elemental and inclusion analysis method, even for non-metallurgists. They will provide an overview of steel cleanliness inspection software that gives full details of your steel sample, including graphic presentation of data in ternary diagrams.
Inclusion Analysis / Steel Cleanliness Determination by Spark OES	

The banner features a collage of scientific images: a circuit board, a glowing orange sphere, a test tube with blue liquid, and a pile of brown grains. The 'FIRST' logo is in large, bold, black letters with white dots. To its right, a blue bar contains the text 'FRONTIERS IN RESEARCH SCIENCE & TECHNOLOGY' in white, and 'newsletter' in white on a dark blue background. The Bruker logo, a blue atom symbol with the word 'BRUKER' in black, is on the right. A small 'ISSUE' label is in the bottom left corner.

**FIRST**

FRONTIERS IN RESEARCH SCIENCE & TECHNOLOGY  
**newsletter**




Like what you learned in this webinar?

Subscribe to Bruker's ***FIRST Newsletter*** to get webinar announcements, fascinating articles, and elemental analysis news delivered right to your inbox!

**Subscribe at:**

<https://www.bruker.com/about-us/register.html>



Innovation with Integrity