

# Elemental Analysis by WDXRF




**Multielement analysis from Beryllium to Uranium in solids and liquids by sequential wavelength dispersive X-ray fluorescence analysis (WDXRF)**




**X-Ray Tube:** Excitation of fluorescence radiation from elements present in the sample

- Low kV for light elements
- High kV for heavy elements



**Primary Beam Filters:** Better detection limits

- Filters out interfering parts of the tube radiation by absorption improving the peak-to-background ratio




**Collimator Masks:** Analysis of different sample sizes

- To cut out radiation from the edge of the sample cup



**Automatic Vacuum Seal:**

- Quick switch of sample types  
Separating sample and goniometer chamber:
- Solid samples: analysis under vacuum
  - Liquid samples: analysis in Helium



**Collimators:** Resolution and sensitivity Collimators form a parallel beam

- Coarse collimators for high sensitivity
- Fine collimators for high resolution



**Analyzer Crystals:** Separation of element specific fluorescence radiation

- Each analyzer crystal covers an element range
- Multilayer optics for enhanced light element analysis



**Detectors:** Detection of the element specific fluorescence radiation

- Light elements are determined with a proportional counter
- Heavy elements are determined with a scintillation counter

