



Product Sheet XRD 37

## VÅNTEC-500 Area Detector for XRD<sup>2</sup> Synchrotron Proven MIKROGAP™ Technology for the Laboratory

The VÅNTEC-500 is a two-dimensional (2-D) detector for X-ray diffraction, based on Bruker AXS' proprietary MIKROGAP™ technology. Due to its huge detector window of 140 mm diameter and its capability of collecting images with both high count rate and excellent resolution, the VÅNTEC-500 is the best detector for all 2-D X-Ray Diffraction (XRD<sup>2</sup>) applications. This includes micro-diffraction, texture analysis, high-throughput diffraction, residual stress determination, Small-Angle X-ray Scattering (SAXS), reciprocal space mapping, etc.

### VÅNTEC-500 - The best detector for XRD<sup>2</sup>

- Huge detection area covering more than 15,000 mm<sup>2</sup>
- Noise-free operation with less than 0.0005 cps/mm<sup>2</sup>
- Variable detector positioning
- Proven radiation-hardness
- Maintenance-free design
- Guaranteed without dead or defective areas



VÅNTEC-500 2-D detector


The VANTEC-500 MIKROGAP detector is the core technology for the next generation of XRD<sup>2</sup> systems by combining the best features of classic X-ray film and electronic single X-ray photon detection. Its huge detector area of more than 15,000 mm<sup>2</sup> grabs as much 2-D information as classic film cameras in just a few seconds. Superb counting statistics is achieved by “seeing” and counting individual photons in real time. The high speed of the VANTEC-500 is ideal for monitoring real time structural and microstructural changes in a material, such as in-situ phase transformation investigation.

Thanks to the proprietary MIKROGAP™ technology, the VANTEC-500 can be operated at very high count rates. Test measurements have shown maximum count rates of 250 Kcps for point-like X-ray signals. At the same time, the VANTEC-500 operates virtually noise-free. The only source of background is cosmic rays and the material’s natural radioactivity, which accounts to less than 5 cps over the

whole detector area. Due to this high local count rate and extremely low background noise, the maximum dynamic range (the ratio of maximum count rate per reflection to background fluctuations in the local area) is excellent: 10<sup>9</sup> times the square root of the collection time in seconds. Hence, the VANTEC-500 is ideal for the analysis of both weakly and strongly scattering samples.


The resistive anode, which is the key feature of the MIKROGAP™ technology, ensures stable detector operation over a wide energy range. The MIKROGAP detector can even tolerate accidental intensive or high energy radiation, as proven by tests at the European Synchrotron Radiation Facility [1]. By design, the VANTEC-500 is also extremely robust with respect to mechanical shock or vibration. Thanks to extensive experience and in-house detector development, the VANTEC-500, like any other detector manufactured by Bruker AXS, is guaranteed to be exemplary and without any defective detector areas.

### MIKROGAP Technology




In a MIKROGAP detector, the grid, the anode, the delay lines and the Xe-CO<sub>2</sub> gas mixture are contained in a sealed vessel with a beryllium entrance window. In the conversion region, the incoming X-ray photons produce primary electrons by gas ionization. These electrons drift towards the anode under a high electric field. The number of electrons is multiplied in the amplification gap. This signal amplification is by far larger than any solid state detector provides thus very weak signals can be detected. The location of each X-ray photon is determined by readout strips through the X and Y delay lines. The key feature of the MIKROGAP technology is the presence of a resistive anode to allow a very thin amplification gap, which greatly improves the local count rate.


### Spatial resolution with MIKROGAP technology



272 µm pixel (512x512)



136 µm pixel (1024x1024)




68 µm pixel (2048x2048)

68 µm pixel size reveals the best spatial resolution, the resolution can be switched by the measurement software

The detector design is optimized with the complete XRD<sup>2</sup> system. Its light weight allows using the detector on both horizontal and vertical diffraction systems. High diffraction angles, e.g. for residual stress measurements, can be achieved thanks to the cone-shaped front side surface. The detector entrance window has a spherical shape to reduce parallax effects, with a radius of the curvature that is designed to match the medium sample-to-detector distance.

In addition, the variable positioning of the VÅNTEC-500 along the track allows for a wide range of sample-to-detector distances to optimize angular coverage and resolution for various applications. In combination with DAVINCI design implemented into the new D8 Diffraction Solution family, the actual distance is even actively recognized by the XRD<sup>2</sup> system and used to automatically pick up the required detector settings.



#### SIZE, the most important feature of 2-D detectors

A large 2-D detector like the VÅNTEC-500 allows the user to:

- Collect several Bragg reflections over a large  $\gamma$  range in one frame
- Measure several polefigures with background correction simultaneously
- Cover broad diffraction peaks of a stressed sample in a single shot
- Enhance the statistics on a spotty diffraction pattern of a micro sample

A large detector window not only enables increased data collection speed, it also provides information that is simply not accessible with a 0-D or 1-D detector (red line). Even small 2-D detectors (orange box) that inherently cover only a limited  $\gamma$  range, will always miss important information.



Technical Specifications	
<b>Sensor Type</b>	Photon-counting, Xe-based MIKROGAP™ detector (US Patent 6,340,819 ) Window size: 140 mm in diameter Number of pixels 2048 x 2048 1024 x 1024 512 x 512 Pixel size 68 µm x 68 µm 136 µm x 136 µm 272 µm x 272 µm
<b>Simultaneous 2θ coverage of a single frame at different sample-to-detector distances</b>	Detector distance 5 cm 10 cm 15 cm 20 cm 25 cm 30 cm 2θ range 83° 56° 42° 33° 27° 23°
<b>Detection Quantum Efficiency</b>	80% for Cu Kα (8.04 keV), residual 20% absorbed by the Be-window
<b>Energy Range</b>	3-15 KeV (Cr, Fe, Co and Cu radiation)
<b>Energy Resolution</b>	20% for Cu Kα (8.04 keV)
<b>Global Counting Rate</b>	Maximum: 1.0 Mcps Linear (10% deviation from linearity): 0.5 Mcps
<b>Local Counting Rate</b>	Maximum per a point-like reflection: 250 Kcps Linear (10% deviation from linearity): 160 Kcps
<b>Background</b>	< 5 cps per whole area (< 0.0005 cps/mm <sup>2</sup> )
<b>Maximum Dynamic Range</b>	10 <sup>9</sup> x √t (t = collection time in seconds)
<b>Radiation Hardness</b>	10 <sup>12</sup> X-ray photons/mm <sup>2</sup> (10 <sup>16</sup> photons in total)
<b>Maintenance</b>	not required
<b>Detector Guarantee</b>	Bruker Detector Guarantee: no defective/dead areas

### References

[1] Nuclear Instruments and Methods in Physics Research A 563 (2006) 172–176

### Bruker AXS GmbH

Karlsruhe · Germany  
Phone +49 721 50997-0  
Fax +49 721 50997-5654  
info.baxs@bruker.com

[www.bruker.com](http://www.bruker.com)