

## SINGLE CRYSTAL X-RAY DIFFRACTION

# APEX5: The central hub of your crystallography programs

Link all your favorite programs for enhanced productivity

### Introduction

APEX5 is the most widely used software suite for single crystal structure determination. It covers the entire process from crystal mounting to CIF report. World-class algorithms ensure reliable and best-in-class results.

Apart from APEX5, there are several other scientific programs developed by the crystallographic community that are specialized in addressing crystallographic problems. These programs, such as ShelXle, PLATON, Mercury, DIAMOND, VESTA, EVAL, JANA2020, or Olex2, offer alternative approaches to crystallography. However, switching between these software packages can be challenging. The poor integration between them requires scientists to repeatedly process the same data and meticulously verify the results for any inconsistencies and errors. These tasks are tedious, error-prone, and redundant, making them time-consuming and unproductive.


### To enhance productivity, it is important to be flexible and agile.

As a flexible and agile crystallographic software suite, APEX5 is designed to meet the requirements of modern scientific research while catering to individual preferences. One of its key strengths is its ability to interface with a range of external crystallographic software packages. The software has been meticulously tested to ensure that data integrity is maintained during the import and export process. This feature significantly boosts productivity, making APEX5 an indispensable tool for scientific research.




## DIAMOND/VESTA

In order to comprehend the arrangement of solid-state compounds, it is essential to conduct a thorough examination of their coordination spheres. Displaying coordination polyhedra serves as a useful method to visualize the coordination environments. DIAMOND and VESTA offer exceptional capabilities in this regard. Simply click on the program link in APEX5 and benefit from these excellent tools for solid-state chemistry.




## JANA2020

Modulated and Quasi-crystals pose an ongoing challenge to most crystallographers. JANA2020 is the gold standard if the periodicity of the crystal under investigation expands to "higher" dimensions. Just click the program link in APEX5 and the multidimensional journey seamlessly continues.



## Olex2

Olex2 and NoSphereA2 can be utilized for building and refining structures. Simply click on the program link icon in APEX5 to continue your work.



## Conclusion


APEX5 serves as the central hub of your crystallographic workflow, providing a direct link to any external software. It allows you to access frequently used programs with a single click, and fully integrates them into the APEX5 interface. This ensures that your latest data is automatically transferred to the third-party program, eliminating the need for manual switching between programs and time-consuming reloading of data. With APEX5 in the center for all your crystallography software needs, you can significantly increase your productivity and streamline your workflow.

## Designed for the Future

APEX5 continues to evolve, driven by invaluable user suggestions and the collective experience of the scientific community. This forward-thinking approach ensures that APEX5 remains at the forefront of cutting-edge crystallography software.

## Empowering the Crystallographic Community

By being a part of the Bruker Users' Group, APEX5 users gain access to a vibrant global family in crystallography. Exchange insights, receive timely responses, and stay informed on the latest developments within the community.



Bruker AXS is continually improving its products and reserves the right to change specifications without notice. Order No. DOC-A86-EX5533 © 2023 Bruker AXS GmbH.

**Bruker AXS**  
info.baxs@bruker.com

bruker.com

**Worldwide offices**  
bruker.com/baxs-offices

**Online information**  
bruker.com/apex5

