



Applications of Elemental Analysis in the Mining and Mineral Resource Industry

Part 2



Welcome to our Minerals, Mining & Geology Webinar Series



Applications of Elemental Analysis in the Mining and Mineral Resource Industry

Part 1: *XRF as a flexible tool for elemental analysis of geological samples*

Dec. 1st 2020: 9:00 AM, 4:00 PM (CET)

→ Get the recording!

Part 2: *Process Monitoring and Grade Control in the Mining and Mineral Resource Industry by XRF*

Dec. 15th 2020: 9:00 AM, 4:00 PM (CET)



Dr. Adrian Fiege

Product Management XRF
Karlsruhe, Germany



Dr. Kai Behrens

Head of XRF Product
Management
Karlsruhe, Germany

Today's Topics



Feedback from Part 1:

- Please tell us more about online-XRF, the S2 KODIAK and typical applications in mining
- Sample preparation as pressed pellet and applications
- Application examples - detection limits and interferences
→ Partially covered



Outline Part 2

- Introduction: When and where to use XRF for Mining operations?
- Sample preparation: Pressed pellets
- Which solution is optimal for my application?
- Application Examples
 - Tungsten and Copper Ore (Online XRF)
 - Iron Ore (EDXRF vs. WDXRF)
 - Nickel Ore (EDXRF vs. benchtop WDXRF)
- Summary & Outlook & Q&A Session



What and when do we analyze in mining operations?



What?

Major, minor and trace elements in all sorts of ores ranging from industrial minerals and commodities, to base metal ore, and to precious metal ores.

-> Focus on decision making elements (grades, tracers, penalty elements)

When?

During exploration

- Individual rock samples
- Drill cores
- Monitoring to survey geological situation, many elements in wide ranges

During mining operation

- Blast holes
- Drill cores
- Quick decisions for operation (no delay), mine planning, major elements

During beneficiation and blending

- Grade control at various steps of the process
- Quick decisions, focus on precision, no further delay of material movement, vital elements

When releasing, shipping, or receiving

- Certify final grade control
- Raw material ID
- Accurate analysis for value determination, highest accuracy



X-ray Fluorescence (XRF) spectrometry

Element range



X-ray Fluorescence (XRF) analysis is qualitative and quantitative method for the determination of element concentrations via excitation of atoms in the sample and detection of the characteristic X-rays.

- High-power WDXRF (4–1 kW):
Be (B) – Am
- Medium WDXRF (400 W):
O (F) – Am
- Modern EDXRF:
C (F) – Am
- Low-power EDXRF:
Na (Mg) – Am

A periodic table of elements where each element's symbol is inside a box. The boxes for elements from Beryllium (Be) to Americium (Am) are highlighted with a yellow bar at the bottom, indicating the range of elements detectable by X-ray Fluorescence (XRF). The highlighted elements include: Be, B, C, N, O, F, Ne, Na, Mg, Al, Si, P, S, Cl, Ar, K, Ca, Sc, Ti, V, Cr, Mn, Fe, Co, Ni, Cu, Zn, Ga, Ge, As, Se, Br, Kr, Rb, Sr, Y, Zr, Nb, Mo, Tc, Ru, Rh, Pd, Ag, Cd, In, Sn, Sb, Te, I, Xe, Cs, Ba, La, Hf, Ta, W, Re, Os, Ir, Pt, Au, Hg, Tl, Pb, Bi, Po, At, Rn, Fr, Ra, Ac, Ce, Pr, Nd, Pm, Sm, Eu, Gd, Tb, Dy, Ho, Er, Tm, Yb, Lu, Th, Pa, U, Np, Pu, Am.


- Element range XRF:
(Be) B to Am
- Concentration range:
Sub-ppm to 100 %

Audience Poll



What are your main goals with elemental analysis – What are you trying to achieve? (Check all that apply.)

- ☐ Check quality of incoming material
- ☐ Increase percentage of recycled material in production
- ☐ Optimize production costs
- ☐ Increase throughput
- ☐ Optimize final product quality
- ☐ Other (specify)



Sample preparation

Pressed pellets vs. Fused beads



Pressed Pellets

Pro's

- Fast & simple sample preparation method
- No loss of volatile elements
- Easily automatable

Con's

- Grain size and matrix effects
- Sample contamination via the grinding vessel



Fused Beads

Pro's

- No Matrix and grain size
- High reproducibility for best accuracy and precision
- Larger calibration range

Con's

- More expensive equipment and consumables
- Higher dilution factor → not suitable for traces

Sample preparation


Pressed pellets



- Fast method
- Easily automatable
- Grain size and matrix effects are not problematic for material-specific calibrations
- But sample contamination via the grinding vessel:
 - **Agate:** SiO_2 99.91%; traces of Al, Na, Fe, K, Mn, Ca, Mg
 - **Corundum ceramic:** Al_2O_3 , traces of K, Na, Si, Ca, Cu, Fe, Mg, Pb, B, Cr, Li, Mn and Ni
 - **Tungsten carbide:** C 6%, Co 6% and W 88%
 - **Chromium steel:** C 1.93%; Cr 13.21%, Traces of Cu, Mn, Mo, Ni, P, Si, S, W

Sample preparation


Pressed pellets



- Grinding was originally carried out with a mortar grinder (0-300 s). The grain size could only be reduced minimally.
- SDM = vibrating disc mill is used today to optimally reduce the grain size
- PP = A press pellet was produced from the SDM sample

Sample preparation

Pressed pellets




Raw meal and limestone

Matrix and grain size effects not compensated

Solution:
Separation of methods

Sample preparation

Pressed pellets




Influence for loose powders and liquids







$$R(\text{kcps}) = f(1/d^2)$$

Δd Anode to sample



Sample not flat:

→ Distance- and shadow-effect







Sample preparation

Grinding



Example: Influence of the grain size



Sample preparation

Thickness of pressed pellet



Z	Element	Line	Energy (keV)	Graphite	SiO2	Fe	Pb
5	B	KA1	0.1833	5.0	1.0	0.3	0.1
9	F	KA1	0.6768	5.0	3.0	0.4	0.3
11	Na	KA1	1.0419	16.0	10.0	0.9	0.9
13	Al	KA1	1.4875	45.0	26.0	2.0	2.0
14	Si	KA1	1.7412	72.0	40.0	4.0	3.0
20	Ca	KA1	3.6910	684.0	88.0	28.0	4.0
56	Ba	LA1	4.4640	989.0	98.0	43.0	5.0
22	Ti	KA1	4.5104	1.3	156.0	47.0	6.0
25	Mn	KA1	5.8981	2.8	338.0	96.0	12.0
26	Fe	KA1	6.4031	3.6	430.0	119.0	15.0
29	Cu	KA1	8.0481	7.1	838.0	28.0	26.0
74	W	LA1	8.3976	8.0	949.0	31.0	29.0
82	Pb	LA1	10.5512	15.0	1.9	57.0	52.0
40	Zr	KA1	15.7749	4.4	6.0	176.0	43.0
42	Mo	KA1	17.4791	5.5	8.0	234.0	47.0
45	Rh	KA1	20.2158	7.3	12.0	355.0	70.0
47	Ag	KA1	22.1630	8.5	16.0	460.0	89.0
56	Ba	KA1	32.1929	13.1	4.0	1.3	242.0

■ = cm

■ = mm

■ = μm

From hand-picked sample to concentrate: Which analytics should I choose?



Hand-held XRF

- Ideal for quick assessment in the field or check of a truck load

Micro-XRF

- Great tool for exploration / drill cores

Online XRF

- Fastest way to check the quality on your conveyor belt and monitor the production process

EDXRF

- Quick and reliable grade control – mainly for major and minor elements

Benchtop WDXRF

- Grade control down to the ppm level

Sequential, full-power WDXRF

- Combining high throughput with accuracy and precision for most elements and concentrations

Multi-channel WDXRF

- Fastest quality check but not flexible regarding changing requirements



S2 KODIAK Online XRF



Online XRF multi-element analyzer for real-time results in mining of base metal ores, and industrial minerals:

- Analysis of major, minor and trace elements
- Real-time analysis of fines and lumps up to > 30 cm diameter
- Covers K – U for lumps
- Covers Si – U for fines
- Handling of dry and humid materials
- Integration time from 1 s upwards
- In mining, mineral beneficiation and blending

S2 KODIAK

Benefits for Mining Applications



- The S2 KODIAK can trigger automated rock sorting with predefined cut-off levels
- Real-time screening of ore grades
- Mineral beneficiation
 - Enhanced recovery rate and higher ore grades
- Blending
 - Produce homogenous ore grades
 - Use highest grade ores sparsely
- Shipping:
 - No penalty claims from ore buyers



S2 KODIAK

Mining Applications



- Exploitation in surface and underground mining
 - Copper
 - Iron
 - Mangan
 - Multimetal
 - Uranium
 - Coal
 - Industrial minerals
- Material sorting and separation of waste rocks when
 - Loading, hauling, crushing, milling, flotation, leaching, concentrating, drying, shipment

S2 KODIAK

Ready for 24/7 operation



- Autonomous operation 24/7
- Integrated UPS, camera and distance sensor
- TCP/IP data transfer to any plant control SW
- Maintenance-free operation, operates with electrical power only
- Optimal occupational health and safety, operates without radioactive sources



- Corrosion resistant, rugged design (V2A)
- Enhanced cooling concept
- Prepared for rough mining environments
- Enclosure protection class
Storage/cleaning mode: IP69K
Acquisition mode: IP65

S2 KODIAK


Installation at tungsten mine



S2 KODIAK, compact, fully radiation safe housing, installed top of the belt
Positioned after first mill, before second mill to analyze feedstock of flotation

Tungsten Ore

Live View of Sample Spectrum



Integration time: 10 s


Power: 3 W

Excitation: 30 kV

Analysis of: K, Ca, Ti, Cr, Mn, Fe, Cu, W,
Zn, Mo

S2 KODIAK

Real versus lab based results



- Excellent match between lab based spot testing and S2 KODIAK
- Better adjustment of flotation parameters based on S2 KODIAK:
 - Drop in concentration of WO_3 on 4.3 and 5.3. is only seen with S2 KODIAK
 - Excesson 6.3. in a short period of time
- **Enhance recovery rates possible**

Project Overview

Application Concentrate



- Analysis of Copper Concentrate material
 - **Sorting of concentrates**
 - High copper content, low contamination content
 - Low copper content, higher contamination content
 - blending process to form homogenized final product
 - The concentrate is typically composed as follows:
 - Cu: 10wt% – 50wt%
 - Fe: 10wt% - 30wt%
 - Contaminants: As, Zn, Bi, Se
 - Material size: Fines


S2 KODIAK

Distance correction

- Material height variation on the belt
- S2 KODIAK operates with
 - Sample height control
 - Integrated ultrasonic sensor
 - Internal Standardization (Argon peak)
 - Software algorithm

Std.Dev.	CuO	Fe2O3
Abs.	0.07 wt%	0.3 wt%
Rel.	0.3 %	1.5 %

Impact of Distance Correction




Dynamic Measurements – Cu Concentrate




- Stability test: one weekend;
ca. 17 h evaluated

	CuO Conc / wt%	Std. Dev. / wt%	Rel. Dev. / %
1.8 cm/s	22.36	0.34	1.6
15 cm/s	22.30	0.34	1.6

Cu Concentrate 1.8 cm/s







Cu Concentrate 15.2 cm/s



Summary

S2 KODIAK Online XRF



-  Within 10 s integration time different Fe grades can be determined and values being used for blending and shipment control
-  Continuous unattended operation possible
-  Arsenic, manganese, calcium possible to integrate w/o additional measurement time
-  Depending on origin different methods can be loaded remotely to the S2 KODIAK



Bruker SIM WDXRF Portfolio


For process control in heavy industries



The S8 LION is our solution for fast, precise and accurate industrial process control in


- Cement
- Industrial minerals
- Base metal mining
- Iron and steel
- Non-Ferrous (Al-Base, Cu-Base)
- Commercial service labs

S8 LION Spectrometer Setup




S8 LION

Rowland Circle Geometry



Spectrum Iron Ore - MEC Dual Mode Data Acquisition



Parallel Measurement of SEC (WDX) and Multi Element Channel
Identification of trace elements, internal backup (second information source)

S8 LION

For central labs and mining service labs



Optimized for highest throughput and optimal cost efficiency

Autonomous operation 24/7

- 90 samples in 6 racks
- 15 samples per rack

In addition to the internal loader

- with 10 positions or,
- with 8 position in combination with online sample feed from automated sample preparation

Add an MEC for extra flexibility



Application Example: Iron Ore




Iron Ore Analysis done to:

- Evaluate the commercial value (mainly Fe total)
- Control the mining process and the refining process
- Highest accuracy and precision in combination with high sample throughput is needed





S8 LION – Iron Ore Pressed Pellets



- Calibration curve for Fe ($R^2 = 0.99983$)
- Elements covered: Fe, Si, Ca, Mn, Al, Ti, Mg, P, S, K, Cr, Co, Ni, Zn

More information with the Multielement Channel MEC



Calibration Details for Fe: Range 34.67 -71.50 %

Fe (Single Element Channel) – LiF 200:

Calibration Std. Dev.: 0.563 %

Fe (Multi Element Channel) – SDD:

Calibration Std. Dev.: 0.503 %

S8 LION – Iron Ore

Repeatability of Pressed Pellets



	Fe [%]	SiO2 [%]	CaO [%]	Mn2O3 [%]	Al2O3 [%]	TiO2 [%]	MgO [%]	P2O5 [%]	SO3 [%]
Rep 01	67.18	1.20	0.045	1.02	0.814	0.033	0.262	0.037	0.102
Rep 02	67.18	1.20	0.045	1.02	0.808	0.033	0.262	0.037	0.102
Rep 03	67.21	1.20	0.045	1.02	0.812	0.033	0.261	0.037	0.105
Rep 04	67.19	1.20	0.045	1.02	0.811	0.034	0.263	0.037	0.106
Rep 05 - 16
Rep 17	67.20	1.20	0.045	1.02	0.812	0.034	0.261	0.037	0.121
Rep 18	67.18	1.21	0.045	1.02	0.815	0.033	0.263	0.037	0.123
Rep 19	67.20	1.20	0.045	1.02	0.811	0.033	0.262	0.037	0.123
Rep 20	67.15	1.20	0.045	1.02	0.812	0.033	0.262	0.037	0.124
Min. [%]	67.14	1.199	0.045	1.018	0.808	0.033	0.260	0.037	0.102
Max. [%]	67.21	1.210	0.045	1.020	0.816	0.034	0.264	0.038	0.124
Mean value	67.17	1.203	0.045	1.019	0.813	0.033	0.262	0.037	0.114
Abs. Std. Dev.	0.022	0.003	0.000	0.001	0.002	0.000	0.001	0.000	0.007
Rel. Std. Dev. [%]	0.033	0.22	0.56	0.058	0.27	1.03	0.28	0.75	6.06

Time per sample: 1:40 min

S8 LION for Incoming Material Inspection Iron Ore/Sinter



- According to ISO 9516 (WDXRF)
- Fused beads
 - At commercial and service labs
 - At central labs
 - At customs
- Maximum 60 s measurement time
- Excellent long-term precision for 200 measurements

Sample	Average (%)	Abs. Std. Dev. (%)	Rel. Std. Dev. (%)
Fe (%)	47.065	0.013	0.028
SiO ₂ (%)	0.673	0.006	0.869
P (%)	0.031	0.001	0.746
Al ₂ O ₃ (%)	4.081	0.014	0.334
Mn (%)	15.348	0.016	0.102
CaO (%)	0.032	0.001	1.867
MgO (%)	0.032	0.002	5.401
TiO ₂ (%)	0.023	0.003	11.705
K ₂ O (%)	0.160	0.001	0.484



S6 JAGUAR

High Performance Benchtop WDXRF



Maintain WDXRF resolution:


- Compact goniometer with high precision gears and closely coupled X-ray beam path
- Optimized analyzer crystals for the entire element range and special applications

Maintain analytical precision and sensitivity:

- Higher power X-ray tube compared to EDXRF (there is no saturation due to single element detection)
- HighSense detection with 2 Mcps count rate
- HighSense XE detector for medium and heavy elements

S6 JAGUAR

Customer feedback





A world-famous manufacturer of refractory and foundry products was looking for a spectrometer for the analysis of main elements in magnesites, etc.,....

- In addition to traces, the analysis of fluorine was also important:
 - Is a large WDXRF spectrometer required?

The S6 JAGUAR showed excellent performance for light and important elements with optimal accuracy, surpassing EDXRF in particular for fluorine.


S6 JAGUAR versus EDXRF F in Foundry Products



Strong overlap of F K α and Fe L α with **EDXRF** leads to medium accuracy and precision: **Min 3.59 % > 3.78 % < Max 4.07 %**


S6 JAGUAR: Optimal resolution, clear separation of both lines, high sensitivity with 400 W power: **Min 3.97 % > 4.03 % < Max 4.07 %**

S6 JAGUAR – Iron Ore Pressed Pellets



- Fe_2O_3 in Iron Ore > 60 % Fe_2O_3

- Fe Ka1 on LiF200 Scan



S6 JAGUAR – Iron Ore

Repeatability of Pressed Pellets




	Al ₂ O ₃ [%]	SiO ₂ [%]	MgO [%]	MnO [%]	Fe ₂ O ₃ [%]	TiO ₂ [%]	P [%]	S [%]
Rep 01	15.91	15.72	0.26	0.02	62.46	0.48	0.118	0.01
Rep 02	15.89	15.70	0.28	0.02	62.68	0.48	0.117	0.01
Rep 03	15.88	15.71	0.25	0.02	62.68	0.48	0.118	0.01
...
Rep 19	...15.88	15.70	0.26	0.02	...62.76	0.48	0.118	0.02
Rep 20	15.85	15.73	0.27	0.02	62.75	0.48	0.116	0.02
Rep 21	15.88	15.72	0.28	0.02	62.61	0.48	0.118	0.02
Mean value	15.89	15.71	0.27	0.02	62.57	0.48	0.12	0.02
Abs. Std. Dev.	0.077	0.032	0.008	0.000	0.379	0.002	0.001	0.004
Rel. Std. Dev. [%]	0.48	0.20	2.95	0.00	0.61	0.46	0.55	19.31

S8 TIGER Series 2 HighSense



- **HighSense Technology** for
 - Ultimate detection limits
 - Optimal sensitivity
 - Best precision
 - Enhanced light element sensitivity
 - Shortest measurement times
- **Optimal instrument configurations** for
 - Industry
 - Metals
 - Minerals & Mining
 - Petrochemistry
 - Cement
 - Chemistry & Pharma
 - Automotive
 - **Academia and Research**
 - Universities
 - Governmental labs
 - Industrial R&D

S8 TIGER Series 2 HighSense



Optimal detection for every element
from Be – Am, ppm to 100%:

- 20 – 60 kV, 5 – 170 mA for 4 kW
- 10 primary beam filters
- 8 collimator masks
- HighSense vacuum seal
- 8 analyzer crystals out of 18 available
- Two detectors with HighSense counting electronics

S8 TIGER Series 2

HighSense: Analyzer Crystals




8 position crystal changer

- More than 18 different crystals available:

optimized **HighSense** crystals:

- XS-CEM: long-term stable, temperature independent
- XS-GE-C: plus 40% intensity for P, plus 20% intensity for S
- XS-PET-C: plus 20 % intensity for Al
- XS-B: plus 100% intensity for B
- XS-N-HS: plus 100% intensity for N
- XS-C: 30% reduced background for C
- XS-100 time optimized measurement
- 25% time savings
- XS-400 plus 35 % more intensity in the range of K to U



S8 TIGER Series 2

HighSense: Counting Electronics



- **HighSense counting electronics**
 - Quick application setup with instant display of PHA scan
 - On-the-fly dead time correction with fast data processing
 - Advanced MCA technology
- Excellent detector linearity
 - FC: up to 4 Mio cps
 - SC: up to 4 Mio cps
 - DynaMatch: up to 13 Mio cps
 - Five times better than conventional WDX spectrometers
- Better analytical precision
- Faster measurements
- Easy application setup due to wide range calibration

GEO-QUANT Iron Ore

Norm compliant Analysis of Iron Ore (ISO 9516)



H																	He
Li	Be											B	C	N	O	F	Ne
Na	Mg											Al	Si	P	S	Cl	Ar
K	Ca	Sc	Ti	V	Cr	Mn	Fe	Co	Ni	Cu	Zn	Ga	Ge	As	Se	Br	Kr
Rb	Sr	Y	Zr	Nb	Mo	Tc	Ru	Rh	Pd	Ag	Cd	In	Sn	Sb	Te	I	Xe
Cs	Ba	La	Hf	Ta	W	Re	Os	Ir	Pt	Au	Hg	Tl	Pb	Bi	Po	At	Rn
Fr	Ra	Ac															

Selection of elements provides

- 21 elements including Na (Na not required in ISO 9516)
- Typically ~16 minutes
- Correct occurring line overlays of the elements
- Matrix effects

Ce	Pr	Nd	Pm	Sm	Eu	Gd	Tb	Dy	Ho	Er	Tm	Yb	Lu
Th	Pa	U	Np	Pu	Am	Sum	Matrix	Compton	Rayleigh				

GEO-QUANT Iron Ore

Concentration Range



	GEO-QUANT Iron Ore		ISO 9516	
	Min [%]	Max [%]	Min [%]	Max [%]
Fe ₂ O ₃	18.9	97.0	54.3	100.0
SiO ₂	0.6	36.3	0.42	13.9
CaO	0.01	15.8	0.03	17.8
Mn ₃ O ₄	0.021	2.0	0.03	1.1
Al ₂ O ₃	0.1	11.5	0.19	6.6
TiO ₂	0.005	10.7	0.03	7.8
MgO	0.02	8.32	0.33	3.3
P ₂ O ₅	0.009	2.7	0.01	1.4
SO ₃	0.007	2.3	0.1	1.5
K ₂ O	0.003	2.6	0.01	0.5
V ₂ O ₅	0.002	0.8	0.003	0.5
SnO ₂	0.0001	0.2	0.001	0.02
Cr ₂ O ₃	0.001	0.2	0.001	0.04
Co ₃ O ₄	0.002	0.2	0.001	0.03
NiO	0.001	0.2	0.014	0.017
CuO	0.008	0.2	0.015	0.08
ZnO	0.002	0.4	0.006	0.21
As ₂ O ₃	0.112	0.1	0.01	0.08
PbO	0.004	0.56	0.02	0.34
BaO	0.0012	0.2	0.04	0.45

- Covers the ISO 9516 application range
- Includes low grade materials
- **QUANT package**
 - 16 certified reference materials (no synthetic standards)
 - Evaluation samples
 - User Manual
 - Preparation Manual



GEO-QUANT Iron Ore

Calibration details – Matrix Correction



Matrix correction with loss eliminated alphas:

- Carbon is set to zero
- New set of alphas has to be calculated
- No further need for LOI determination



Elemental Analysis of Ni-Laterite Ore


Fast and reliable with the S2 PUMA Series 2



Nickel Laterite (low grade nickel ore) will likely be the **dominant source of nickel** in the foreseeable future.

EDXRF is used for grade control in mining operations and smelters. It allows, e.g.,

- accurate determination of Ni
- analysis of waste rock elements
 - Ca, Mg, Fe, Mn
- evaluation of traces
 - Ti, Cr, Co, Cu, Zn
- Bruker is a trusted EDXRF supplier with many S2 PUMA installations at Ni-Laterite mines
- Time-to-result, ease-of-use, and flexibility are major benefits of the S2 PUMA with XY Autochanger



Elemental Analysis of Ni-Laterite Ore

Fast and reliable with the S2 PUMA Series 2




The New S2 PUMA Series 2

Enhanced Speed

- 3 x Higher count-rate, thanks to HighSense technology
- 2 x Higher pump-rate
- Faster Software (up to 90%)

Enhanced Usability & Flexibility

- Improved user interface
- New features, such as Advanced/Basic mode and a Drift Monitoring tool
- ... and much more
- Setting thresholds for "High" and "Low" grade ore allows to quickly evaluate the results

A screenshot of the S2 PUMA Series 2 software interface, titled "SPECTRA-RESULTS MANAGER". It displays a table of analysis results for Nickel Grade and Ni (%) across multiple sample IDs. The interface includes a menu bar, a toolbar with icons for various functions, and a status bar at the bottom.

Sample ID	Nickel Grade	Ni (%)
UNKNOWN4	LOW	1.95
UNKNOWN4	LOW	1.96
UNKNOWN4	LOW	1.96
UNKNOWN4	LOW	1.96
194	HIGH	2.10
194	HIGH	2.11
194	HIGH	2.10
194	HIGH	2.11

Elemental Analysis of Ni-Laterite Ore

Fast and reliable with the S2 PUMA Series 2





- The new hardware features of the S2 PUMA Series 2 allow to further optimize the analytical settings for Ni-laterite analysis, resulting in better performance and shorter measurement times
- The improved workflow, the faster software, and the higher pump rate decrease the time-to-result even more
- **Analytical procedure:**
 - 25 kV, no beam filter, fixed current, 60 s counting time
 - Pressed pellet preparation



Elemental Analysis of Ni-Laterite Ore

Fast and reliable with the S2 PUMA Series 2



- Calibration curve for Ni in Ni-Laterite (left) and Ni KA1 peaks (right)
- $R^2 = 0.99887$

Elemental Analysis of Ni-Laterite Ore

Fast and reliable with the S2 PUMA Series 2



Sample ID	Time-to-Result [min:sec]	Ni (%)	MgO (%)	Al ₂ O ₃ (%)	SiO ₂ (%)	CaO (%)	TiO ₂ (%)	Cr ₂ O ₃ (%)	MnO (%)	Fe ₂ O ₃ (%)	Co (PPM)	Cu (PPM)	Zn (PPM)
Rep-1	01:37	1.388	24.05	2.094	48.24	0.32	0.03	0.77	0.21	14.02	375	32	119
Rep-2	01:39	1.393	24.19	2.077	48.17	0.32	0.03	0.77	0.21	14.00	377	31	125
Rep-3	01:39	1.403	24.00	2.102	48.10	0.32	0.03	0.77	0.22	14.05	416	35	119
Rep-4	01:31	1.383	24.08	2.102	48.19	0.32	0.03	0.77	0.21	14.00	376	32	116
Rep-5	01:39	1.393	24.00	2.107	48.13	0.32	0.03	0.77	0.22	14.05	374	33	116
Rep-6	01:32	1.382	24.28	2.086	48.35	0.32	0.03	0.77	0.22	13.96	400	34	124
Rep-7	01:39	1.39	24.18	2.102	48.26	0.32	0.03	0.77	0.22	13.97	369	33	122
Rep-8	01:32	1.384	24.09	2.094	48.19	0.32	0.02	0.77	0.21	14.03	336	33	123
...
Rep-120	01:31	1.407	23.84	2.088	47.85	0.33	0.03	0.78	0.22	14.11	357	33	120
Average measured value	01:31	1.391	24.09	2.101	48.19	0.32	0.03	0.77	0.22	14.02	358	33	120
Abs. standard deviation		0.007	0.08	0.012	0.10	0.00	0.00	0.00	0.00	0.04	21	1	5
Rel. Standard deviation		0.50	0.32	0.57	0.21	1.46	8.00	0.26	2.25	0.29	5.77	3.31	4.03

- **Excellent analytical performance in ~1.5 min** (time-to-result)
(most EDXRF spectrometers require 4-5 min)

XRF – The right tool for all tasks

S8 LION:

- Simultaneous instrument with shortest measurement times (<2 min), 4 kW, highest precision

S8 TIGER

- Sequential WDXRF for high demanding grade control. Analytical flexibility for major, minors and traces, optimal settings for each elements

S6 JAGUAR

- Benchtop WDXRF: lower throughput but similar performance compared to full-size WDXRF, better spectral resolution than ED

S2 PUMA

- Versatile tool for minor and major elements for quick decisions in mobile and small labs

S2 KODIAK

- On-belt analysis for quickest automated decisions along the production of high grade materials



Q&A




Any questions?



Thanks for your time and interest!





Innovation with Integrity